

Product datasheet for **MR209552L4V**

Ddx46 (BC026492) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ddx46 (BC026492) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ddx46
Synonyms:	MGC31579
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC026492
ORF Size:	1878 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209552).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC026492.1
RefSeq Size:	2545 bp
RefSeq ORF:	1880 bp
Locus ID:	212880
Cytogenetics:	13 B1
Gene Summary:	Plays an essential role in splicing, either prior to, or during A complex formation. [UniProtKB/Swiss-Prot Function]



[View online »](#)