

Product datasheet for **MR209558L4V**

Lcp1 (NM_008879) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lcp1 (NM_008879) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lcp1
Synonyms:	AW536232; D14Erttd310e; LCP-1; Pls2; pp65
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008879
ORF Size:	1884 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209558).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008879.3
RefSeq Size:	3748 bp
RefSeq ORF:	1884 bp
Locus ID:	18826
UniProt ID:	Q61233
Cytogenetics:	14 39.63 cM
Gene Summary:	Actin-binding protein. Plays a role in the activation of T-cells in response to costimulation through TCR/CD3 and CD2 or CD28. Modulates the cell surface expression of IL2RA/CD25 and CD69.[UniProtKB/Swiss-Prot Function]



[View online »](#)