

Product datasheet for **MR209568**

Smpd1 (NM_011421) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Smpd1 (NM_011421) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Smpd1
Synonyms:	A-SMase; ASM; aSMase; Zn-SMase
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR209568 representing NM_011421
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCCCACCACAGAGCATCATCCGGCCAGGACCCTCAGAGCCGGCTGGGAGCAGAGACTGGAGAGGT
 CCTTACCGGCACCCAGAGTGGGACTCCTTTGGATGGGGCTGGGCTTGGCGCTGGTTCTGGCTCTGTTTGA
 CTCACCGTTCTTTGGGTTCTGCCAGAGCTTATCCTCTTCTTGAAGGCCATTCTGTCAAATTCAGT
 GCCATAGCGCCGCGCTCCAGAGTGCCTTTGGGTGGCAGAACCTCACTTGCCCGCCTGCAAAGTCTTAT
 TCACTGCTCTCAACCATGGGCTGAAGAAGGAGCCCAATGTGGCACGGGTAGGCTCTGTGGCCATCAAGAT
 ATGCAAGATGCTGAACATAGCACCCTAGATGTGTGCCAGTCAAGCCGTCATCTCTTTGAGGACGATGTG
 GTGGAGGTGTGGACACGTTCACTTCTGAGCCCATCAGAGGCTTGTGGCTTGTCTTGGGCTCCTCTGTG
 GACTGGGACATCTTTTCGACTTGAACATCTCTTTGCCATCAGTGCCGAAGCCACCCCAAAGCCACC
 GAGCCACCGCCCGGTGCCCTGTAGCCGTGCTCTTCTTACTGATCTACTGGGACCATGAG
 TACCTGGAGGGCAGACCCCTTACTGTGAGATCCACTTTGCTGCCGCGGGGGTCCGGATGGCCACCCA
 ACTCCCAGAAGGGGGCTGGGTTCTGGGGCGAGTACAGCAAGTGTGACCTGCCCTTGCACGCTAGAAA
 CCTGTTGAAAGGACTGGGCCCTGCCGCCCTTTTGAATGGTGTACTGGACGGGAGATATCCCTGCCCAT
 GATGTCTGGCAACAGTCTCGCAAGATCAGCTGAGGGCCCTGACCACCATCACAGACCTCGTGAGGAAGT
 TCTTGGGCCAGTGCCAGTGTACCCTGCTGTGGCAACCATGAGAGCACTCCTGTCAATGGCTTCCCTCC
 TCCCTTATAAAGGGAACCAATCTTACAATGGCTCTATGAAGCCATGGCAAGGATGGGAACCTGG
 TTACCAGCTGATGCCCTCACACCCTAAGAATGGGGCTTCTATGCCCTACCCACGCCCTGGCCTCC
 GCCTCATCTCTCAATATGAATTTTGTCCCGTGAGAATTTTGGCTTTGATCAACTCCACAGATCC
 TGCTGGACAACCTCAGTGGCTGGTGGAGGAGTTCAGGCTGCTGAGAATCGAGGAGACAAAGTGCATATA
 ATTGGCCACATCCCTCCAGGACATTGTCTTAAGAGCTGGAGCTGGAATTATTACAAAATCATAGCCAGT
 ATGAAAACACTCTGGCCGGTCACTTCTTTGGCCACACTCACGTGGATGAGTTTGGATCTTCTATGATGA
 GGAAACTCTGAGCCGACCACTAGCTGTAGCCTTCTGGCGCCAGTGTACAACCTTTATCAACCTTAAC
 CCTGGCTACCGAGTTTACCAAATAGATGGAACTACCCCGAAGCTCTCATGTGGTCTGGACCACGAGA
 CCTACATCTTGAATCTGACCCAGGCCAACGCAGCAGGAGGCACACCAAGCTGGAAGCGCCTACAGGGC
 TCGAGAAACCTATGGACTTCCAGATGCGATGCCTGCCTCTGGCACAACCTGGTCTACCGCATGAGGGAT
 GATGAGCAACTCTCCAGACCTTCTGGTTCTCTACCATAAGGGCCACCCACCTCAGAGCCCTGCGGCA
 CACCCTGCCGCCCTGGCCACTCTGTGTGCCAGCTCTCAGCCCAGCAGACAGCCCTGCTTTGTGTGCCA
 CTTGATGCCCAATGGGAGCTCCAGATGCTAATAGTTGTGGTACGCCCCCTGTTGTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR209568 representing NM_011421
 Red=Cloning site Green=Tags(s)

MPHHRASSGQDHLRAGWEQRLERSLPAPRVGLLWMGLGLALVLALFDSTVLWVPARAYPLPSEGHVSKFS
 AIAPPLQSAFGWQNLTPACKVLFALNHGLKKEPNVARVGSVAIKICKMLNIAPLDVQCQSAVHLFEDDV
 VEVWTRSVLSPSEACGLLLGSSCGHWDIFSTWNI SLPSVPKPPKPPSPPAPGAPVSRVFLDLDLHWDHE
 YLEGTDPYCADPLCCRRGSGWPPNSQKAGAGFWGEYSKCDLPLRRTLESLLKGLGPAGPFEMVYWTGDIPAH
 DVWQSRQDQLRALTTITDLVRKFLGPVPVYPAVGNHSTPVNGFPPPF IKGNQSSQWL YEAMAKAWEPW
 LPADALHTLRIGGFYALTPRPLRLISLNMNFC SRENFWLLINSTDPAGQLQWLVEELQAAENRGDKVHI
 IGHIPPGHCLKSWSWNYKIIARYENTLAGQFFGHTHVDFEIFYDEETLSRPLAVAFAPSATTFINLN
 PGYRVYQIDGNYPGSSHVLDHETYILNLTQANAAGGTPSWKRLYRARETYGLPDAMPASWHNLVYRMRD
 DEQLFQTFWFLYHKGHPPEPCGTPCRLATLCAQLSARADSPALCRHLMPNGSLPDANRLWSRPLLC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9014_d08.zip
 Restriction Sites: SgfI-MluI
 Cloning Scheme:



ACCN: NM_011421

ORF Size: 1881 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_011421.2](#), [NP_035551.1](#)

RefSeq Size: 2430 bp

RefSeq ORF: 1884 bp

Locus ID: 20597
UniProt ID: [Q04519](#)
Cytogenetics: 7 55.9 cM
MW: 70.4 kDa
Gene Summary: Converts sphingomyelin to ceramide. Also has phospholipase C activities toward 1,2-diacylglycerolphosphocholine and 1,2-diacylglycerolphosphoglycerol.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209568