

Product datasheet for **MR209624**

Hspa2 (NM_001002012) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hspa2 (NM_001002012) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Hspa2
Synonyms:	70kDa; Hsp70-2; HSP70.2; HSP70A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR209624 representing NM_001002012
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGC
 GCC

ATGTCTGCCCGCGGCCGGCTATCGGCATCGACCTGGGCACCACTTACTCGTGCCTGGGGTATTCCAAC
 ATGGCAAGGTGGAGATCATCGCCAACGACCAGGTAACCGCACCAACCCAGCTACGTGGCCTTCACTGA
 CACCGAGCGTCTCATCGGCGACGCCAAGAACCAAGTGGCCATGAACCCACAAAACCATCTTCGAC
 GCCAAGCGGCTGATCGGACGGAAGTTCGAGGACGCCACAGTGCAGTCCGATATGAAGCACTGGCCGTTC
 GAGTGGTGAAGGTTGGGAAGCCAAAGTGCAGGTGGAATAAAAGGGGAGATGAAGACCTTCTTCCC
 TGAGGAGATTTCTCCATGGTCTCACTAAGATGAAGGAGATCGCGAAGCCTACCTGGGGGCAAGGTG
 CAGAGCGCAGTCATCACTGTTCTGCCTATTTCAACGACTCGCAGCGCCAGGCCACCAAGGATGCGGGCA
 CCATCACCGCCCTCAACGTGTTGCGCATCATCAATGAGCCACAGCAGCGGCCATCGCCTACGGCCTGGA
 TAAGAAGGGCTGTGCGGGCGGCGAGAAGAAGTGTCTCATCTTTGACCTGGGCGGGGCGACCTTCGATGTG
 TCCATCTGACCATCGAGGATGGCATCTTTGAGGTGAAGTCCACGCGCCGCGCATACCCACCTGGGTGGCG
 AAGACTTCGACAACCGTATGGTCAAGCCCTGGCGGAGGAGTTCAAGCGCAAACAAGAAGGACATTGG
 GCCCAACAAGCGCGCTGTGCGCCGGTGCACCCGCTGCGAGCGCGCTAAGCGCACCTGAGCTCGTCC
 ACGCAGGCCAGCATAGAGATCGACTCGCTCTACGAGGGCGTGGATTTCTACACGTCCATCACCCGCGCC
 GCTTCGAGGAGCTCAACGCCGATCTCTCCGAGGGACCTGGAGCCGGTGGAAAAGGCGCTGCGCGATGC
 CAAGCTAGACAAGGGCCAGATCCAGGAGATAGTGTGTTGGGCGGCTCAACCCGCATCCCTAAGTCCAG
 AAGCTCCTGCAAGATTTCTTCAACGGCAAGGAGCTGAACAAGAGCATTAAATCCCGACGAGGCGGTGGCCT
 ACGGCGCCGCTGTGACGCGGCTATCTCATCGGCGACAAGTCGAAAAATGTCAGGATCTGCTGTTACT
 CGACGTGACTCCATTGTCGCTCGGCATCGAAACAGCTGGCGGTGTCATGACCCCACTCATCAAGAGAAAC
 ACCACGATCCCCACCAAGCAGACGACGACCTTCACTACCTACTCAGACAACCAGAGCAGCGTGTGGTGC
 AAGTGTACGAGGGCGAACGGGCCATGACCAAGGACAATAACCTCTTGGGCAAGTTCGACCTGACTGGGAT
 CCCCCAGCACCCGTTGGGTCCCCAGATCGAGGTACCTTTGACATCGATGCCAACGGCATCCTTAAC
 GTCACTGCTGCCGACAAGAGCACCGGTAAGAAAATAAAATCACCATACCAACGACAAGGGTGGGCTGA
 GCAAAGACGACATTGACCGGATGGTGCAGGAGGCGGAGCGGTACAATCGGAAGATGAAGCAAATCGCGA
 TCGCGTGGCAGCAAAAATGCGGTGGAGTCTATACCTACAACATCAAGCAGACCGTGAAGACGAGAAA
 CTGAGGGCAAGATTAGCGAGCAGGACAAAACAAGATCCTCGACAAGTGTGAGGAGGTGATCAACTGGC
 TTGACCGAAACCAGATGGCAGAGAAAGATGAGTACGAACAAGCAGAAAAGAGCTTGAGAGAGTGTGCAA
 CCCCATCATCAGCAAACTTTACCAAGGCGGTCCAGGCGGCGGCGCTCCTCTGGAGGGCCACCATCGAG
 GAAGTGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >MR209624 representing NM_001002012
 Red=Cloning site Green=Tags(s)


```
MSARGPAIGIDLGTTYSCVGVFQHGKVEIANDQGNRTTPSYVAFDTERLIGDAAKNQVAMNPTNTIFD
AKRLIGRKFEDATVQSDMKHWPFRVVSSEGGKPKVQVEYKGMKTFPPEEISSMVLTKMKEIAEAYLGGKV
QSAVITVPAYFNDSQRQATK DAGTITGLNVLRIINEPTAAAIAYGLDKKGCAGGEKNVIFDLGGGTFDV
SILTIEDGIFEVKSTAGDTHLGGEDFDNRMVSHLAEFEFKRKHKKDIGPNKRAVRRRLTACERAKRTLSSS
TQASIEIDSLYEGVDFYTSITRARFEELNADLFRGTLEPVEKALRDAKLDKGQIQEIVLVGGSTRIPKIQ
KLLQDFFNKELNKSINPDEAVAYGAAVQAAILIGDKSENVQDLLLDDVTPLSLGIETAGGVMTPLIKRN
TTIPTKQTQFTTYS DNQSSVLVQVYEGERAMTKDNNLLGKFDLTGIPPAPRGV PQIEVTFDIDANGILN
VTAADKSTGKENK IITNDKGR LSKDDIDRMVQEAERYKSEDEANRDRVA AKNVAVSYTYNIKQTV EDEK
LRGKISEQDNKILDKCQEVINWLD R NQMAEKDEYEHKQKELERVCNPIIISKLYQGGPGGGGSSGGPTIE
EVD
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

AscI-MluI

Cloning Scheme:



ACCN: NM_001002012

ORF Size: 1899 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001002012.1 , NP_001002012.1
RefSeq Size:	2550 bp
RefSeq ORF:	1902 bp
Locus ID:	15512
UniProt ID:	P17156
Cytogenetics:	12 33.73 cM
MW:	69.6 kDa
Gene Summary:	<p>Molecular chaperone implicated in a wide variety of cellular processes, including protection of the proteome from stress, folding and transport of newly synthesized polypeptides, activation of proteolysis of misfolded proteins and the formation and dissociation of protein complexes. Plays a pivotal role in the protein quality control system, ensuring the correct folding of proteins, the re-folding of misfolded proteins and controlling the targeting of proteins for subsequent degradation. This is achieved through cycles of ATP binding, ATP hydrolysis and ADP release, mediated by co-chaperones. The affinity for polypeptides is regulated by its nucleotide bound state. In the ATP-bound form, it has a low affinity for substrate proteins. However, upon hydrolysis of the ATP to ADP, it undergoes a conformational change that increases its affinity for substrate proteins. It goes through repeated cycles of ATP hydrolysis and nucleotide exchange, which permits cycles of substrate binding and release (By similarity). Plays a role in spermatogenesis (PubMed:24557841). In association with SHCBP1L may participate in the maintenance of spindle integrity during meiosis in male germ cells (PubMed:24557841).[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR209624