

## Product datasheet for **MR209631L4V**

### **Tmem201 (NM\_001025106) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	Tmem201 (NM_001025106) Mouse Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	Tmem201
<b>Synonyms:</b>	9430091K04; AV028368; D4Ertd429e; Net5; RP23-38317.1; Samp1
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-mGFP-P2A-Puro (PS100093)
<b>Tag:</b>	mGFP
<b>ACCN:</b>	NM_001025106
<b>ORF Size:</b>	1905 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(MR209631).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_001025106.1</a> , <a href="#">NP_001020277.1</a>
<b>RefSeq Size:</b>	4620 bp
<b>RefSeq ORF:</b>	1905 bp
<b>Locus ID:</b>	230917
<b>Cytogenetics:</b>	4 80.15 cM
<b>Gene Summary:</b>	Involved in nuclear movement during fibroblast polarization and migration (PubMed:22349700). May recruit Ran GTPase to the nuclear periphery (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)