

Product datasheet for **MR209646**

Pip5k1c (NM_001146687) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pip5k1c (NM_001146687) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pip5k1c
Synonyms:	A1115456; A1835305; Pip5klgamma
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR209646 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGCTAGAGGTGCCGACGAGGCGGAGAGCGCCAAGCGGGGCGGTGACGGCGGAAGCGCCCTGGT
CTGCGGAGAGTGGGGCGGCGCAGGTATGACCCAGAAGAAGGCTGGCCTTGCAGAGGCCCTTGGTGAC
AGGGCAGCCAGGCCCTGGCCATGGGAAGAAGCTGGGCCACCGAGGCGTGGATGCGTCGGGAGAGACTACG
TATAAGAAGACCACCTCGTCCACCCTGAAGGGCGCCATCCAGCTGGGGATCGGGTACACGGTGGCAACC
TGAGCTCTAAGCCGAGCGGGACGTGCTCATGCAGGACTTCTACGTGGTGGAGAGCATCTTCTCCCCAG
TGAAGGGAGCAACCTCACCCCTGCCACCATTCCAGGATTTCCGCTTCAAGACCTATGCACCTGTTGCC
TTCCGCTACTTCCGGGAGCTCTTGGCATCCGTCCAGATGACTATTTGACTCCCTGTGCAATGAGCCAC
TCATCGAGCTCTCAACCCCGGGGCCAGCGGCTGTCTTCTACGTACCAGCGACGACGAGTTCATCAT
CAAGACTGTCATGCACAAGGAGGACAGATTCTGCAGAAGCTGCTGCCTGGCTACTACATGAATCTCAAC
CAGAACCACGGACGCTGCTGCCAAGTTCTATGGGCTGTACTGCGTGCAGTCTGGTGGCAAGAATCC
GCGTGGTGGTCAATGAACAATGTGCTGCCCGTGTGCTCAAAATGCACCTTAAGTTCGACCTCAAGGGCTC
CACGTACAAGCGCAGGGCCAGCAAGAAGGAGAAGGAGAGCCTGCCACCTACAAGGACCTGGACTTC
ATGCAGGACATGCCCGAGGGGCTGCTGCTGGACTCCGACACCTTTGGCGCCCTGGTCAAGACGCTGCAGC
GAGACTGCCTGGTGTGGAGAGCTTCAAGATAATGGACTACAGCCTGCTGCTGGGCGTGCACAACATCGA
TCAGCAGGAGCGAGAGCGCCAGGCCGAGGGCGCCAGAGCAAGGCGGATGAGAAGCGGCCGTGGCCAG
AAGGCTCTGTATCCACGGCCATGGAGTCTATCCAGGGCGGAGCTGCCCGTGGGGAGGCCATTGAGACAG
ATGACACGATGGGTGGGATTCCAGCAGTGAATGGGCGCGGGGAGCGACTGCTCCTCACATCGGGATCAT
TGATATTCTGCAGTCTACAGTTCATCAAGAAGTTAGAACACACCTGGAAGGCCCTCGTCCATGATGGG
GACACTGTCTCAGTCCACCGGCCAGTTCATGCAGAGCGCTTCTTCAAGTTCATGAGCAGCACGGTGT
TCCGGAAGAGTTCCTCCCTGAAGTCTCTCCATCCAAGAAAGGACGTGGTGCCTGCTGGCGGTCAAACC
CCTGGGGCCCACTGCTGCCTTCTCAGCCAGCCAGATCCCAGCGAGAGAGAGGATGTGCAGTATGACCTG
CGGGGGCCCGCAGTACCCACGCTTGAGGATGAAGGCCGGCCTGACCTCCTGCCCTGCACCCACCGT
CCTTTGAGGAAGCCACCACAGCCTCCATCGCCACCACCCTGTCGTCACCTCCCTCTCCATCCCAGAGCG
GTCCCCTTCAGATACATCGGAGCAGCCCGGTACAGGCGGCGCACGCAGTCTTACGGCCAGGATGGCCGG
CCCCAGGAGGAGCCCATGCGGAAGACCTGCAGAAGATAACTGTGCAGGTGGAGCCAGTGTGCGGTGTGG
GGTTGTCCCAAGGAGGAGGGTGCAGGAGTGGAGGTCCCCCATGTGGGGCATCGGCTGCAGCCTCTGT
GAAATAGACGCTGCCAGCCAGGCTCAGAGCCTGCCAGCCAGGCTCAGATGAGGAGGATGCACCCCTCT
ACAGACATCTATTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR209646 protein sequence
Red=Cloning site Green=Tags(s)

```
MELEVPDEAESAEAGAVTAEAAWSAESGAAAGMTQKKAGLAEAPLVTGQPGPGHGKLLGHRGVDASGETT
YKKTTSSTLKGAIQLGIGYTVGNLSSKPERDVLMDQDFYVVESIFFPSEGSNLTPAHHFQDFRFTYAPVA
FRYFRELFGIRPDDYLYSLCNEPLIELSNPGASGSVFYVTSDDDEFIIKTMHKEAEFLQKLLPGYMNLN
QNPRTLTPKFYGLYCVQSGGKNIRVVVMNNVLPVVVMHLKFDLKGSTYKRRASKKEKEKSLPTYKDLDF
MQDMPEGLLLDSDTFGALVKTLQRDCLVLESFKIMDYSLLLGVHNIDQQERERQAEGAQSKADEKRPVAQ
KALYSTAMESIQGGAARGEAIETDDTMGGIPAVNGRGERLLLHIGIIDILQSYRFIKKLEHTWKALVHDG
DTVSVHRPSFYAERFFKMSSTVFRKSSSLKSSPSKKGREGALLAVKPLGPTAAFSASQIPSEREDVQYDL
RGARSYPTLEDEGRDLLPCTPPSFEEATTASIATTLSSLSIPERSPSDTSEQPRYRRRTQSSGQDGR
PQEEPHAEDLQKITVQVEPVCVGVVPEEGAGVEVPPCGASAAAASVEIDAASQASEPASQASDEEDAPS
TDIYF
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001146687

ORF Size: 1908 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146687.2](#), [NP_001140159.1](#)

RefSeq Size: 4222 bp

RefSeq ORF: 1908 bp

Locus ID: 18717


UniProt ID: [O70161](#)

Cytogenetics: 10 C1

MW: 69.5 kDa

Gene Summary: Catalyzes the phosphorylation of phosphatidylinositol 4-phosphate (PtdIns4P) to form phosphatidylinositol 4,5-bisphosphate (PtdIns(4,5)P2). PtdIns(4,5)P2 is involved in a variety of cellular processes and is the substrate to form phosphatidylinositol 3,4,5-trisphosphate (PtdIns(3,4,5)P3), another second messenger. The majority of PtdIns(4,5)P2 is thought to occur via type I phosphatidylinositol 4-phosphate 5-kinases given the abundance of PtdIns4P. Participates in a variety of cellular processes such as vesicle mediated transport, cell adhesion, cell polarization and cell migration. Together with PIP5K1A is required for phagocytosis, but they regulate different types of actin remodeling at sequential steps. Promotes particle attachment by generating the pool of PtdIns(4,5)P2 that induces controlled actin depolymerization to facilitate Fc-gamma-R clustering. Mediates RAC1-dependent reorganization of actin filaments. Required for synaptic vesicle transport. Controls the plasma membrane pool of PtdIns(4,5)P2 implicated in synaptic vesicle endocytosis and exocytosis. Plays a role in endocytosis mediated by clathrin and AP-2 (adaptor protein complex 2). Required for clathrin-coated pits assembly at the synapse. Participates in cell junction assembly. Modulates adherens junctions formation by facilitating CDH1 trafficking. Required for focal adhesion dynamics. Modulates the targeting of talins (TLN1 and TLN2) to the plasma membrane and their efficient assembly into focal adhesions. Regulates the interaction between talins (TLN1 and TLN2) and beta-integrins. Required for uropodium formation and retraction of the cell rear during directed migration. Has a role in growth factor- stimulated directional cell migration and adhesion. Required for talin assembly into nascent adhesions forming at the leading edge toward the direction of the growth factor. Negative regulator of T-cell activation and adhesion. Negatively regulates integrin alpha-L/beta-2 (LFA-1) polarization and adhesion induced by T-cell receptor. Together with PIP5K1A has a role during embryogenesis and together with PIP5K1B may have a role immediately after birth. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209646