

Product datasheet for **MR209658**

Gclc (NM_010295) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gclc (NM_010295) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gclc
Synonyms:	D9Wsu168e; Ggcs-hs; GLCL-H; Gclc
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR209658 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGGCTGCTGTCCAAGGCTCGCCACTGAGCTGGGAAGAGACCCAGCGCCACGCCACCAGTGC
 GGCACGGCATCCTCCAGTTCCTGCAGATCTACCACGCAGTCAAGGACCGGCACAAGGACGTGCTCAAGTG
 GGGTGACGAGGTGGAGTACATGTTGGTGTCTTCGATCATGAAAAAGGAAAGTCCAGTTGTTACTGAAT
 GGAGGCGATGTTCTTGAGACTCTGCAGGAGAAGGGGGAGAGGACAAACCCCAACCATCCGACCCTCTGGA
 GACCAGAGTATGGGAGTTACATGATCGAAGGAACGCCTGGCCAGCCCTACGGAGGAACGATGTCTGAGTT
 CAACACTGTGGAGGCAATATGAGGAAACGCCGGAAGGAGGCTACTTCTGTACTAGGAGAGCATCAGGCT
 CTCTGCACCATCACTTCATCCCCAGACTAGGCTGTCTGGATTCACACTGCCAGAACACAGACCCCAACC
 CAGAGGAAGGAGGCGCATCAAAGTCCCTCTTCCAGATGAAGCCATAAACAAGCACCCCGCTTCGG
 TACTCTAAACAAGAAACATCCGGCATCGGAGAGGAGAAAAGTTGTCATCAATGTGCCAATATTCAAGGAC
 AAGAATACACCATCTCCATTTGTAGAAACATTTCCAGAAGATGCGGAGGCATCAAAGGCTTCTCAGCCAG
 ACCATATCTACATGGATGCCATGGGATTCGGGATGGGCAACTGCTGTCTCCAGGTGACATTTCAAGCCTG
 CAGCATATCTGAGGCAAGATACCTTTATGATCAGCTGGCCACTATCTGCCAATTGTTATGGCTTTGAGT
 GCTGCATCTCCATTTTACCGAGGCTACGTGTCAGACATTGATTGTGCTGGGGAGTGATTCTGCCTCTG
 TAGATGATAGAACACGGGAGGAGAGAGGCTGGAGCCTCTGAAGAACAATCGCTTTAGGATCAGTAAATC
 TCGGATGACTCAATAGATAGCTACCTGTCTAAGTGTGGGAGAAAGTACAATGACATCGACCTGACCATC
 GATAAGGAGATCTACGAGCAGCTCCTGGAGGAAGGCATCGATCATCTTCTGGCACAGCACGTTGCTCATC
 TCTTTATTAGAGACCCACTCACCTTTTTGAAGAGAAAAATACACCTGGATGATGCCAACGAGTCTGACCA
 TTTTGAGAAATTCAGTCTACAAACTGGCAGACAATGAGGTTTAAAGCCTCCTCCTCCAAACTCAGATATT
 GGATGGAGAGTAGAGTTCGACCAATGGAGGTGCAGTTAACAGACTTTGAGAACTCTGCCTATGTGGTAT
 TCGTGGTACTGCTCACTAGGGTATCCTCTCATACAAAGTACTTCTCATTCCGCTGTCCAAGGTTGA
 CGAGAACATGAAAGTGGCCAGAAAGCGAGATGCTGTCTTGACGGGATGTTTTATTTAGGAAAGACATC
 TGCAAAGGCGGCAACGCTGTGGTGGACGGGTGCAGCAAGGCCAGAGCAGCTCTGAGCCAGCTGCAGAGG
 AGTACACCCTGATGAGCATAGACACCATCATCAATGGGAAGGAAGGGTATTTCTGGACTCATCCCCAT
 TCTGAACTCCTATCTGGAGAACATGGAAGTTGATGTGGACACCCGATGCAGTATTCTGAACTACCTAAAG
 CTAATTAAGAAGAGAGCATCTGGAGAATAATGACTGTTGCCAGGTGGATGAGAGAGTTTCATTGCAAACC
 ATCCTGACTACAAGCAAGACAGTGAATAACTGATGAAATCAACTATAGCCTCATTTGGAAGTGCAACCA
 AATTGCAGATGAGCTGTGTGAATGTCCAGAGTTACTTGGATCAGGATTTAGAAAAGCGAAGTACAGTGGA
 GGTAAAAGCGACCCTTCAGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR209658 protein sequence
 Red=Cloning site Green=Tags(s)


MGLLSQGSPLSWEETQRHADHVRRHGILQFLQIYHAVKDRHKDVLKWGDEVEYMLVSFDHENRKVQLLN
 GGDVLETLQEKGERTNPNHPTLWRPEYGSYMIIEGTPGQPYGGTMEFNTVEANMRKRKEATSVLGEHQA
 LCTITSFPRLGCPGFTLPEHRPNPEEGGASKSLFFPDEAINKHPRFGTLTRNIRHRRGEKVVINVPFKD
 KNTPPSFVETFPEDAASKASQPDHIYMDAMGFGMGNCCLQVTFQACSISEARYLYDQLATICPIVMALS
 AASPFYRGYVSDIDCRWGVISASVDDRTREERGLEPLKNNRFRIKSKRYDSIDSYLSKCGEKYNDIDLTI
 DKEIYEQLLEEGIDHLLAQHVAHLFIRDPLTLFEEKIHLDDANESDHFENIQSTNWQTMRFKPPPPNSDI
 GWRVEFRPMEVQLTDFENSAYVVFVLLTRVILSYKLDLFIPLSKVDENMKVAQKRDAVLQGMFYFRKDI
 CKGGNAVVDGCSKAQSSSEPAEEYTLMSIDTIINGKEGVFPGLIPILNSYLENMEVDVDRCSILNYLK
 LIKKRASGELMTVARWMREFIANHPDYKQDSVITDEINYSLIWKCNIADDELCEPELLGSGFRKAKYSG
 GKSDPSA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_010295

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_010295.1](#), [NM_010295.2](#), [NP_034425.1](#)

RefSeq Size: 3442 bp
RefSeq ORF: 1914 bp
Locus ID: 14629
UniProt ID: [P97494](#)
Cytogenetics: 9 43.36 cM
MW: 72.6 kDa