

Product datasheet for **MR209665L4V**

Cstf2t (BC026995) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cstf2t (BC026995) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cstf2t
Synonyms:	tCstF-64, 64kDa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC026995
ORF Size:	1911 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209665).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC026995 , AAH26995
RefSeq Size:	3583 bp
RefSeq ORF:	1913 bp
Locus ID:	83410
Cytogenetics:	19 C1
Gene Summary:	May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs.[UniProtKB/Swiss-Prot Function]


[View online »](#)