

Product datasheet for **MR209666**

Limk2 (NM_001034030) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Limk2 (NM_001034030) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Limk2
Synonyms:	Limk-2; whe
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>MR209666 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGCGCTGGCGGGTATGAAGCCTGGAGGTGTCGGGGCTGTGGCACCTATGTTCCCTGAGCCAGA
GGCTCTATAGGACTGCCAATGAAGCCTGGCATGGCTCCTGCTCCGGTGTCCGAATGCCAGGAATCCCT
CACCAACTGGTACTATGAGAAGGATGGGAAGCTGTACTGCCACAAGGACTACTGGCCAAGTTTGGAGAA
TTCTGCCATGGGTGCTCCCTGCTGATGACGGGGCCAGCCATGGTGGCTGGGGAGTTTAAGTACCACCCAG
AGTGCTTTGCCTGCATGAGTTGCAAGGTGATCATTGAGGATGGAGATGCATATGCCCTGGTGCAGCACGC
TACTACTACTGTGGAAAGTGCCACAACGAGGTGGTGTGGCACCCATGTTTAAAAGACTGTCTACAGAG
TCTGTCCAGGATCAGCTGCCCTACTCTGTACGCTCATCTCCATGCCTGTACCACCGAGTGTAGACCGG
GCTTCTCGGTGACTGTGAAAAGTGCCTCTTCTAACTATGCCACCACCGTTCAAGTAAAAGAGGTCAACCG
GATGCACATCAGTCCCAACAACCGAAATGCCATCCACCCTGGGGACCGCATCTTGGAGATCAATGGGACC
CCTGTCCGTAATCTGCGAGTAGAGGAGGTGGAGGATGCAATAAAGCAGACAAGCCAGACACTTCAGCTGC
TGATTGAACATGACCCTGTGCCCCAGCGCTGGACCAGCTACGGCTAGATGCCCGACTTCCTCCCCACAT
GCAGAGACCCGGACACACGCTCATGCTCAGCACCTTGGACACCAAGGAGAATCAGGAGGGGACACTGAGG
AGACGTTCTCTGAGGCGCAGTAACAGCATCTCCAAGTCTCCGGGCCCCAGCTCCCCAAGGAGCCCTGC
TCCTCAGCCGGGACATAAGCCGCTCAGAATCCCTCCGCTGCTCTAGCAGTACTCCCAGCAGATCTTCCG
GCCCTGCGACCTGATCCACGGGGAGTCTGGGAAGGGCTTCTTTGGGCAGGCCATCAAGGTGACTCAC
AAAGCCACAGGCAAAGTGATGGTCATGAAGGAGTTGATCCGCTGTGATGAGGAAACCCAGAAGACTTCC
TGACTGAGGTGAAAGTGATGCGGAGCCTGGACCACCTAATGTGCTCAAGTTCATCGGAGTCTGTACAA
GGACAAGAAGCTGAATCTGCTGACCGAGTACATCGAGGGAGGCACACTCAAAGACTTCCTGCGCAGTGTG
GACCCGTTCCCTGGCAACAGAAGTCAAGTTTGCCTAAAGGCATCTCCTCTGGAATGGCCTATTTACT
CCATGTGCATCATCCACCGGACCTGAACCTCACACAAGTGTCTCATCAAATGGACAAGACAGTGGTAGT
AGCCGACTTTGGGCTGTACGGCTTATAGTAGAAGAGAGGAAAAGGCCTCCAGTAGAGAAGGCCACCACC
AAGAAACGCACCTTACGCAAGAGTGACCGCAAAAAGCGCTACTGTGGTGGGAAACCCCTACTGGATGG
CCCCGAGATGCTGAATGGCAAGAGCTACGATGAGACGGTGGATGTCTTCTTTTGGGATCGTTCTCTG
TGAGATCATTGGGCAGGTATATGCGGATCCTGATTGCCTGCCCGCACACTGGACTTTGGCCTCAATGTT
AAGCTTTTCTGGGAGAAGTTTGTCCCAACAGACTGTCCCCAGCCTTCTCCCCCTGGCTGCTATCTGCT
GCAAACCTAGAGCCTGAAAGCAGACCGGCATTCTCAAAGCTGGAGGACTCGTTCGAGGCGCTCCCTGTT
CTTGGGGGAGCTGGCCATCCCCTGCCAGCAGAGCTGGAAGACCTGGACCACACTGTGAGCATGGAGTAT
GGCCTAACCCGGGACTCGCCACCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR209666 protein sequence
 Red=Cloning site Green=Tags(s)


MAALAGDEAWRCRCGCTVYVPLSQRLYRTANEAWHGSCFRCSECQESLTNYYEKDGKLYCHKDYWAKFGE
 FCHGCSLLMTGPAMVAGEFKYHPECFACMSCKVIIEDGDAYALVQHATLYCGKCHNEVVLAPMFERLSTE
 SVQDQLPYSVTLISMPATTECRRGFSVTVESASSNYATTVQVKEVNRMHISPNRNNAIHPGDRILEINGT
 PVRTL RVEEVEDA IKQTSQTLQLLIEHDPVQRLDQLRLDARLPPHMQSTGHTLMLSTLDTKENQEGTLR
 RRSLRRSNSISKSPGPSSPKPELLLSRDISRSESLRCSYSQQIFRPCDLIHGEVLGKGFQAIKVTH
 KATGKVMVMKELIRCDEETQKFTLEVKVMRSLDHPNVLKFIGVLKDKKLNLLTEYIEGGTLKDFLRV
 DPFWQQKVRFAKGISSGMAYLHSMCIHRDLNSHNCLIKLTKTVVADFGLSRLIVEERKRPPVEKATT
 KKRTLKSDRKKRYTVVGNPYWMAPEMLNGKSYDETVDVFSFGIVLCEIIGQVYADPDCLPRTLDFGLNV
 KLFWEKFVPTDCPPAFFPLAAICCKLEPESRPAFSKLEDSFEALSLFLGELAIPLPAELEDLHTVSM EY
 GLTRDSPP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001034030

ORF Size: 1917 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001034030.1](#), [NM_001034030.2](#), [NP_001029202.1](#)

RefSeq Size:	3898 bp
RefSeq ORF:	1356 bp
Locus ID:	16886
UniProt ID:	<u>O54785</u>
Cytogenetics:	11 A1
MW:	72.2 kDa
Gene Summary:	Displays serine/threonine-specific phosphorylation of myelin basic protein and histone (MBP) in vitro.[UniProtKB/Swiss-Prot Function]