

Product datasheet for MR209689L4

Fubp1 (BC023780) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fubp1 (BC023780) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Fubp1
Synonyms:	9530027K12Rik; D3Ertd330e; FBP; Fubp; Fubp4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209689).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 51886

Cytogenetics: 3 77.05 cM

Gene Summary: Regulates MYC expression by binding to a single-stranded far-upstream element (FUSE) upstream of the MYC promoter. May act both as activator and repressor of transcription (By similarity).[UniProtKB/Swiss-Prot Function]