

## Product datasheet for **MR209689L4V**

### Fubp1 (BC023780) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fubp1 (BC023780) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fubp1
Synonyms:	9530027K12Rik; D3Ertd330e; FBP; Fubp; Fubp4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC023780
ORF Size:	1920 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209689).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC023780.1</a>
RefSeq Size:	2425 bp
RefSeq ORF:	1922 bp
Locus ID:	51886
Cytogenetics:	3 77.05 cM
Gene Summary:	Regulates MYC expression by binding to a single-stranded far-upstream element (FUSE) upstream of the MYC promoter. May act both as activator and repressor of transcription (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)