

Product datasheet for **MR209716L4V**

Tm9sf4 (NM_133847) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tm9sf4 (NM_133847) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tm9sf4
Synonyms:	AA986553; AU045326; B930079E06; mKIAA0255
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133847
ORF Size:	1929 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209716).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133847.3 , NP_598608.2
RefSeq Size:	3895 bp
RefSeq ORF:	1932 bp
Locus ID:	99237
UniProt ID:	Q8BH24
Cytogenetics:	2 H1
Gene Summary:	Associates with proteins harboring glycine-rich transmembrane domains and ensures their efficient localization to the cell surface.[UniProtKB/Swiss-Prot Function]



[View online »](#)