

## Product datasheet for **MR209738**

### **Srrm3 (NM\_021403) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Srrm3 (NM_021403) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Srrm3
Synonyms:	2900083I11Rik; AB041661; SRm300-like; Srrm2l
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR209738 representing NM\_021403  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGC**C

ATGTCCTCCACCGTGAACAACGGGGCTACCGGCATGCCAGCTCCACCGGATGCTGCTAATGGCTTCCCGC  
AGCCCGCGCCTCGTCCGGTCTTGGCCGCGCGGGAGGAGAGCTGCGCGCCGCGAGCCGGCTTGGT  
GAAGCGCGCGCACCGGGAGATCCTGGACCAGAACGCAAGCGGCGCTGGAGCTCAAGTGCATGGAGTTG  
CAGGAGATGATGGAGGAGCAGGGTACTCAGAGGAGGAGATCCGCGAGAAGGTCGGGACCTTCCGGCAGA  
TGCTGATGGAGAAGGAAGCGTACTGACCCGGGAGGACCGGCCTGGGGCCACATCGTGGCAGAGACCC  
GCGACGAATGGAGGGCTTGGAGCCGGCCTGGAGTACCCACCCTTCGATGAGGACGACGGCCGGTGGAT  
TGCGACTGCCCGTCTCCTGCTACCGCGGGCACCGGGGTACAGGACCAAGCACTGGTCCAGCAGCTCTG  
CATCGCCCCCTCCAAGAAGAAGAAGAAAAGAAAGGCAGTCATCGGAGGAGCCGAAAAAGAGGAGACT  
AGAGTCAGAATGCAGCTGTGGGAGTGTCCCCCTGCGCAAGAAGAAGAATGTGAAGAAACACCGT  
CGAGATCGGTCTGATTCTGGGTCCCGGAGGAAGAGACGATACAGATCGAGAAGCCTCAAGAGCAAACGAA  
AAGAAAAAACAAGAGCGAAAGAGGCCTCACACTGAGTCTCCAGGCCGAAGTTCCACCACCACAGCAG  
CGCCAGCTCCCACAGCCCTCAATGTCTCACACTACAGTATTCTGGATCACCCAGCAGACTGAGCCCC  
AAGCATCGAGATGACGGGCGGAAGACAGGCAGCCAGCATCCAGCGGGAGCCGGTCCGATCCCCATCGG  
GCGGCAGCGGATGGGGTTCGCCCCAGCAGAACGGAGGCAGCAGGCAGGCAAGCGGAGCGCACGGGGTGC  
CCCCGGCTCGGCGCACAGCCCGCCGATAAGCCAGCTCGCCCCGGGCTGCGACAAGCGGGCCGCGCGG  
CCCACGCCTCCCGCGCTGGGAAGGACAGCCAGAGTCCCCGCTCCGCGCCGTCGTCGACGGGCCGCGGG  
GCCGCGCGGGGAGGGACAGCCAGGCGCCGGTCCGGCAGCGGAGGAGAGACGATCGCGTCTCAGGC  
CAACGCGCCCCGCGCAGGGGACGCCGGGTACAAAGCCCGCGCCGCTCGGGGCTCGTCGCGCTCGCTC  
AGCGGGGCACACTCCAGCAGCGACTCCGGTGGCGGCGCCCGGGCCCGGGCCCGAGCCCTGCTCCGAGC  
GAGGTCACGGCGGACACGGGAAACGGGCAAGGAGCGCCCCCGCGCGCACGCCCGGCTAGCACCTCGCC  
ATCGCCTGGCACGCGGGCCGACGGGGCGGCCAGAGGGCAACAGCTCATCGCGCAGCCGGTCCCCAT  
CAGGGCTCCTGGAGCTCCAGCCGCTCGCCCTCAAATCCACTCGCGCTCCCAGACAAAAGGACCCGCA  
GCCCCAGCCTCTACCGTCGCCAAGAAACCTCTTGGCAGGGACAAGGACAGCGAAGGCCGTGCGAGGCA  
CGCCGAGGCTGAGGCCGCCGACCCGGGTCGCTCCCGCAGCTACTCGCCATCCGAAAGCGGCGCGG  
GACTCGCCGAGTTTCATGGAGCCGAGACGCATCACAGTCCCGCAAGCGTCCCATCCCCTACTACCGGC  
CCAGCCCCTCCTTCTCCAGCTGCTTGGACCCGACTACTCAAGCCGGAGCCACAGCCGGAGCCCCAG  
CCCTGGCCACAGCCATGGAAGCTACAGCAGCCGAGTACCGGACTCGCAGCCGCTCCTGCAGCGCTCC  
AGGAGCCGAGTCTAGCTACCACAGCCGGAGCAGCTCTGAGAGTGGGGCTTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR209738 representing NM\_021403  
 Red=Cloning site Green=Tags(s)


MSSTVNNGATGMPAPPDAANGFPQPGASSGSWPRAEEELRAAEPGLVKRAHREILDHERKRRVELKCMEL  
 QEMMEEQGYSEEEIRQKVGTFRQMLMEKEGVL TREDRPGAHIVAE TPRRMEGLEPGL EYPPFDEDDGPVD  
 CDCPVS CYRHRGYRTKHWSSSASPPPKKKKKKKGSHRRSRKKRRLESECSCGSASPLRKKKKNVKKHR  
 RDRSDSGSRRKRRYRSRSLKSKRKEKNERKRPHTESPGRRFHHHSSASSHSPSMSSHYSDSGSPSRLSP  
 KHRDDGKRKTGSQRSSGSRSPSPSGSGWGSPQQNGGSRQRSGAHGGRPGSAHSPDPKSSPRACDKAAAA  
 PTPPARGKDSQSPRSAPSSQGRGGAAGGTARRRRRRRRRRRSRSANAPRRRGRRTKPAPPRGSSRSL  
 SGAHSSSDSGGGAPGPGPEPCSERGHGGH GKRAKERPPRAR PASTSPSPGTRRRGGPEGNSSSRSPGPH  
 QGSWSSSRSPSKSHSRSPDKRTRSPSLSPSPKKPLGRDKDSEGRARHAEAE AARTRRRSYSPIRKR  
 RRRDPSPFMEPRRITSARKRPIPYRPS PSSSSSCLSTDYSSRSHSRSPSPGHSHGSYSSRSHGTRSRSCSAS  
 RSRSPSYHSRSESSEGGF

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9108\\_d05.zip](https://cdn.origene.com/chromatograms/mm9108_d05.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_021403

**ORF Size:** 1944 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_021403.3](#), [NP\\_067378.2](#)

**RefSeq Size:** 3970 bp

**RefSeq ORF:** 1947 bp

**Locus ID:** 58212


**UniProt ID:** [Q80WV7](#)

**Cytogenetics:** 5 G2

**MW:** 71.1 kDa

**Gene Summary:** May play a role in regulating breast cancer cell invasiveness. May be involved in RYBP-mediated breast cancer progression.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR209738