

Product datasheet for **MR209742L3V**

Slc23a2 (NM_018824) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc23a2 (NM_018824) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc23a2
Synonyms:	AI844736; mKIAA0238; NBTL1; Slc23a1; SVCT2; YSPL2; YSPL3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018824
ORF Size:	1944 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209742).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018824.2 , NP_061294.2
RefSeq Size:	6487 bp
RefSeq ORF:	1947 bp
Locus ID:	54338
UniProt ID:	Q9EPR4
Cytogenetics:	2 F2
Gene Summary:	Sodium/ascorbate cotransporter. Mediates electrogenic uptake of vitamin C, with a stoichiometry of 2 Na(+) for each ascorbate (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)