

Product datasheet for **MR209762**

Ddx17 (NM_001040187) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ddx17 (NM_001040187) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ddx17
Synonyms:	2610007K22Rik; A430025E01Rik; AI047725; C80929; Gm926; p7; p72
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>MR209762 representing NM_001040187
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCGTGGAGCGGCTTTGGGGATCGGGATCGTGACAGGGACCGTGGAGGTTTGGAGCAAGAGGTGTA
GTGGGCTTCCCCCTAAGAAGTTTGGTAATCCTGGGGAGCGGTTACGAAAAAGAAGTGGATTTGAGTGA
ACTCCCTAAATTTGAGAAGAATTTTATGTTGAGCATCCAGAAGTAGCAAGACTGACTCCGTATGAGTT
GATGAGCTACGGCGTAAGAAAGAGATTACAGTGAGAGGGGAGATGTTTGTCCAAAACCTGTCTTTGCC
TCCATCATGCTAACTTTCCACAGTATGTGATGGATGTGCTGATGGATCAGCACTTTACAGAACCCACTCC
CATTCACTGCCAGGGATTTCTTTGGCTCTTAGTGGCAGGGATATGGTTGGCATTGCACAGACTGGCTCT
GGGAAGACATTGGCGTATTTGCTGCCTGCGATTGTTTATATAAACACCAGCCATACTTGGAAAGAGGAG
ATGGTCCAATTTGTCTAGTGTGCTGCTACCAGAGAGCTTGCAGCAGGTGCAGCAGGTGGCTGACGA
TTATGGAAAATGCTCCAGGTTGAAGAGTACGTGCATTTACGGAGGTGCTCCTAAAGTCCCAAATTCGA
GACTTGGAAAGAGGTGTTGAGATTTGCATAGCCACTCCTGGGCGCCTAATAGATTTCTGGAGTCAGGAA
AGACAAACCTTCGCCGATGACTTACCTTGTGTTGGATGAGGCTGACCGGATGCTTGATATGGGCTTTGA
GCCCCAGATCCGTAATTTGTTGATCAAATCAGGCCTGACCGGCAGACACTGATGTGGAGTGCAACCTGG
CCGAAGGAAGTGAGGCAGCTTGCAGAGGATTTCTGCGGGACTACACCCAGATCAATGTGGGCAATCTGG
AGCTGAGTGCCAACCACAACATCCTACAGATTGTAGATGTCTGTATGGAGAGTAAAAAGACCACAAATT
GATCCAGCTGATGGAGGAAATCATGGCGGAAAAGGAAAATAAGACTATAATATTTGTGGAGACAAAGAGG
CGCTGTGATGACCTCACACGAAGAATGCGCAGAGATGGTTGGCCTGCTATGTGTATCCATGGAGACAAGA
GTCAACCAGAAAGAGATTGGGTACTTAATGAGTCCGATCTGGAAGGCTCCTATCCTCATTTGCCACGGA
TGTAGCCTCCCGTGGGCTAGATGTGGAAGATGTCAAGTTTGTATCAACTACGATTATCCAAACAGCTCA
GAGGATTATGTTACCGTATTGGCCGAACGCGCCGAGCACCAACAAGGGCACTGCCTATACTTTCTTTA
CCCCGGGCAACCTGAAGCAGGCTAGAGAGCTGATCAAAGTATTGGAAGAGGCAATCAAGCCATCAATCC
AAAATTGATGCAGCTTGTGGACCACAGAGGTGGCGGGAGGAGGGGGAGGCCGCTCACGATACCGGACT
ACTTCTTCAGCCAACAATCCCAATCTGATGTATCAGGACGAGTGTGACCGGAGGCTTCGAGGGGTCAAGG
ATGGTGGCCGGAGAGATTCTACAAGCTACAGGGATCGTAGTGAAACCGATAGAGCCAGTTATGCTAATGG
CAGTGGCTATGGAAGCCAAATTTCTGCCTTTGGGGCACAAGCAGGCAATACACCTATGCTCAAGGCACC
TATGGGGCAGCTGCCTATGGCACCAGTGGCTACACGGCGCAGGAGTATGCTGCTGGCACTTACGGGGCGA
GCAGCACTGCCTCAGCAGGGAGGAGCTCTCAGAGCTCCAGCCAGCAGTTTGTGGGATAGCCGATCTGG
GCAGCAGCCACAGCCACTGATGTACAGCAGTTTGCACAGCCTCCAGGAGCTACCAATATGATAGGCTAC
ATGGGGCAGACTGCTTACCAGTACCCTCCCCCTCCCCCTCCCCCTCCTCCATCTCGCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR209762 representing NM_001040187
 Red=Cloning site Green=Tags(s)


MRGGGFGDRDRDRDRGGFGARGGSGLPPKFGNPGERLRKKKWDLSELPKFEKNFYVEHPEVARLTPYEV
 DELRRKKEITVRGGDVC PKPVFAFHANFPQYVMDVLMQHFTEPTPIQCQGFPLALSGRDMVGIAQTGS
 GKTLAYLLPAIVHINHQPYLERGDGPICLV LAPTRELAQQVQQVADDYGKCSRLKSTCIYGGAPKGPQIR
 DLERGVEIC IATPGRLIDFLESKTNLRRCTYLVLDEADRMLDMGFEPQIRKIVDQIRPDRQTLMWSATW
 PKEVRQLAEDFLRDYTQINVGNELELSANHNILQIVDVCMESEKDHKLIQLMEEIMAEKENKTIIFVETKR
 RCDDLTRMRRDGWPAMCIHGDKSQPERDWLNEFRSGKAPIL IATDVASRGLDVEDVKFVINYDYPNSS
 EDYVHRIGRTARSTNKGTA YTFFTPGNLKQARELIKVLEEANQAINPKLMQLVDHRGGGGGGGSRRYRT
 TSSANNPNLMYQDECDRRLRGVKDGGRRDSTSYRDRSETDRASYANGSGYSPNSAFGAQAGQYTYAQGT
 YGAAAYGTSGYTAQEYAAGTYGASSTASAGRSSQSSSQFSGIGRSGQPQLMSQQFAQPPGATNMIGY
 MGQTAYQYPPPPPPPPSRK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:
 Cloning Scheme:

SgfI-MluI

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

Plasmid Map:


ACCN: NM_001040187

ORF Size: 1950 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001040187.1](#), [NP_001035277.1](#)

RefSeq Size: 4766 bp

RefSeq ORF: 1953 bp

Locus ID: 67040

UniProt ID: [Q501J6](#)

Cytogenetics: 15 E1

MW: 72.8 kDa

Gene Summary: This gene encodes the mouse homolog of human DEAD box polypeptide 17. DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD). RNA helicases of the DEAD-box family are involved in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Alternative splicing of this gene results in several transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]