

## Product datasheet for **MR209785L4V**

### **Pdzrn3 (BC010329) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	Pdzrn3 (BC010329) Mouse Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	Pdzrn3
<b>Synonyms:</b>	1110020C07Rik; AI429718; AL023082; LNX3; SEMACAP3; Semcap3
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-mGFP-P2A-Puro (PS100093)
<b>Tag:</b>	mGFP
<b>ACCN:</b>	BC010329
<b>ORF Size:</b>	1962 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(MR209785).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">BC010329</a> , <a href="#">AAH10329</a>
<b>RefSeq Size:</b>	2921 bp
<b>RefSeq ORF:</b>	1964 bp
<b>Locus ID:</b>	55983
<b>Cytogenetics:</b>	6 D3
<b>Gene Summary:</b>	E3 ubiquitin-protein ligase. Plays an important role in regulating the surface level of MUSK on myotubes. Mediates the ubiquitination of MUSK, promoting its endocytosis and lysosomal degradation. Might contribute to terminal myogenic differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)