

## Product datasheet for **MR209802L2V**

### **Mx2 (NM\_013606) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Mx2 (NM_013606) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mx2
Synonyms:	AI528743; Mx-2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_013606
ORF Size:	1968 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209802).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013606.1</a> , <a href="#">NP_038634.1</a>
RefSeq Size:	2424 bp
RefSeq ORF:	1968 bp
Locus ID:	17858
UniProt ID:	<a href="#">Q9WVP9</a>
Cytogenetics:	16 57.51 cM



[View online »](#)

**Gene Summary:**

This gene encodes a member of the Mx protein family of large GTPases, and functions in the innate immunity system. Interferon alpha/beta treatment or viral infection induces expression of this protein, which subsequently accumulates in the cytoplasm and inhibits viral replication. It has been shown to confer resistance to vesicular stomatitis virus, a member of the rhabdovirus family, but not to influenza virus. This gene produces a functional protein in some feral mouse strains, whereas some inbred mouse strains including the strain of the reference genome, C57BL/6J, contain a frameshift that inactivates this gene product. [provided by RefSeq, Jul 2008]