

Product datasheet for MR209811L3

Acadvl (NM_017366) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Acadvl (NM_017366) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Acadvl
Synonyms:	VL; vlcad
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209811).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_017366
ORF Size:	1971 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017366.2
RefSeq Size:	2205 bp
RefSeq ORF:	1971 bp
Locus ID:	11370
UniProt ID:	P50544
Cytogenetics:	11 42.96 cM
Gene Summary:	This gene encodes a homodimeric mitochondrial flavoprotein and is a member of the acyl-CoA dehydrogenase family. Members of this family catalyze the first step of fatty acid beta-oxidation, forming a C2-C3 trans-double bond in a FAD-dependent reaction. As beta-oxidation cycles through its four steps, each member of the acyl-CoA dehydrogenase family works at an optimum fatty acid chain-length. This enzyme has its optimum length between C16- and C20- acylCoA and localizes to the inner mitochondrial membrane (unlike related acyl-CoA dehydrogenases). In mice, deficiency of this gene can cause ventricular arrhythmias as well as fasting and cold intolerance. [provided by RefSeq, Nov 2012]

Product images:



Circular map for MR209811L3