

## Product datasheet for **MR209813LV**

### Ripk1 (NM\_009068) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ripk1 (NM_009068) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ripk1
Synonyms:	D330015H01Rik; Rinp; RIP; RIP-1; Rip1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009068
ORF Size:	1971 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209813).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009068.3</a> , <a href="#">NP_033094.3</a>
RefSeq Size:	4512 bp
RefSeq ORF:	1971 bp
Locus ID:	19766
UniProt ID:	<a href="#">Q60855</a>
Cytogenetics:	13 14.01 cM



[View online »](#)

**Gene Summary:**

Serine-threonine kinase which transduces inflammatory and cell-death signals (programmed necrosis) following death receptors ligation, activation of pathogen recognition receptors (PRRs), and DNA damage (PubMed:12654725, PubMed:19590578). Upon activation of TNFR1 by the TNF-alpha family cytokines, TRADD and TRAF2 are recruited to the receptor (By similarity). Phosphorylates DAB2IP at 'Ser-728' in a TNF-alpha-dependent manner, and thereby activates the MAP3K5-JNK apoptotic cascade (By similarity). Ubiquitination by TRAF2 via 'Lys-63'-link chains acts as a critical enhancer of communication with downstream signal transducers in the mitogen-activated protein kinase pathway and the NF-kappa-B pathway, which in turn mediate downstream events including the activation of genes encoding inflammatory molecules (By similarity). Polyubiquitinated protein binds to IKBKG/NEMO, the regulatory subunit of the IKK complex, a critical event for NF-kappa-B activation (By similarity). Interaction with other cellular RHIM-containing adapters initiates gene activation and cell death (By similarity). RIPK1 and RIPK3 association, in particular, forms a necrosis-inducing complex (By similarity).[UniProtKB/Swiss-Prot Function]