

Product datasheet for **MR209819L3V**

Abcg2 (NM_011920) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Abcg2 (NM_011920) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abcg2
Synonyms:	ABC15; ABCP; AI428558; BCRP; Bcrp1; MXR; MXR1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011920
ORF Size:	1974 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209819).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011920.3 , NP_036050.1
RefSeq Size:	2493 bp
RefSeq ORF:	1974 bp
Locus ID:	26357
UniProt ID:	Q7TMS5
Cytogenetics:	6 27.82 cM



[View online »](#)

Gene Summary:

The membrane-associated protein encoded by this gene is included in the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the White subfamily. Alternatively referred to as a breast cancer resistance protein, the human protein functions as a xenobiotic transporter which may play a major role in multi-drug resistance. This protein likely serves as a cellular defense mechanism in response to mitoxantrone and anthracycline exposure. [provided by RefSeq, Jul 2008]