

Product datasheet for **MR209826L3V**

Gpsm1 (BC071197) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gpsm1 (BC071197) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gpsm1
Synonyms:	1810037C22Rik; Ags3; AW107933
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC071197
ORF Size:	1974 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209826).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC071197 , AAH71197
RefSeq Size:	3474 bp
RefSeq ORF:	1976 bp
Locus ID:	67839
Cytogenetics:	2 A3



[View online »](#)

Gene Summary:

Guanine nucleotide dissociation inhibitor (GDI) which functions as a receptor-independent activator of heterotrimeric G-protein signaling. Keeps G(i/o) alpha subunit in its GDP-bound form thus uncoupling heterotrimeric G-proteins signaling from G protein-coupled receptors. Controls spindle orientation and asymmetric cell fate of cerebral cortical progenitors. May also be involved in macroautophagy in intestinal cells. May play a role in drug addiction. [UniProtKB/Swiss-Prot Function]