

Product datasheet for **MR209830**

Kif22 (NM_145588) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kif22 (NM_145588) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Kif22
Synonyms:	AU021460; C81217; Kid; Kif22a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR209830 representing NM_145588
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAGCTTAAGGGCTAAGACGTGTCCGACGACGGGAGATGGCGTCGGCCACCTCAGGACCTGGCCGCT
 GTGTAAGCAAAGGAGGGCTTGGCCGGCGTCCGCCCTTAGCTCGAGTAAGGGTGGCTGTGCGGCTGCGCCC
 GTTTATGGATGGAGAGACAGAAGCGAAGGAGCTCCCTGTGTCCGAGCCATAGACAGCTGCTCTCTTGAA
 GTGGCTAACTGGAAGAAATACCAGGAGACGCTCAAATATCAGTTTGATGCCTTCTATGGCGAGAAGAGCA
 CTCAGCAGGAAGTCTATGTAGTTTCAGTACAGCCCATCTAAGGCACTTGTGGAAGGGCAGAATGCCAG
 TGTACTTGCTATGGGCCTACTGGGCAGGGAAGACACACAATGCTGGCAGCCAGAGCAACCTGGA
 GTGATTCCTCGGGCTCTCATGGACCTCTGCAACTAGCAAGGGAAGAGAGTGTGAGGGACGGCCATGGG
 ACGTTTCTGTTGCTATGCCTATTTAGAGATCTACCAGGAAAAGTATTAGACCTCTTGATCCTGCATC
 AGGAGACCTCGTGATCCGCGAAGACTGCCGAGGGAACATCCTGATCCAGGCCTCACACAGAAGCCCATC
 ACCAGTCTCTGACTTCGAGCAGCACTTCCTTCCAGCCAGTAGAAATCGAGCTGTAGGAGCCACCCGGC
 TTAACCAGCGCTCTTCCCGTAGTCACGCAGTGTCTTGGTCAAGGTAGATCAGCGTGAACGTTTGACTCC
 ATTTCCGACGCGGGAAGGAAACTCTACCTTATTGATTTGGTGGTTCAGAGGACAACCGTCGCACAGGC
 AACCAGGGCATTCCGGCTCAAAGAGAGCGGAGCCATCAACACCTCCCTCTTTGACTGGGCAAAGTGGTGG
 ATGCATTAACCAGGGCCTCCCTCGCATACCATACCGGACAGCAAGTCACTCGTCTGCTGCAGGACTC
 TCTGGGAGGCTCAGCTCATAGCATCCTCATTGCCAACATTGCTCCTGAAAGACGTTTTTACCAGGATACA
 ATTTCAGCATTGAACTTCACTGCTAGGTCCAAGGAGGTGATTAACCGGCCTTCCACCAATGAGAGTCTAC
 AGCCTCATGCCTTGGCACCTGTTAAGCTGTCTCAGAAAGAACTGCTAGGCCCGTCAGAGGCAAAGAAGGC
 CAAAGGCCCTGAAGAAGAATCAACTGGGAGCCCTGAGTCTACAGCAGCTCCAGCGTCTGCCTCCAGAAA
 CTAAGCCTCTTACAGAAGCTAAGCAATATGGATCCAGCCATGTTGGAGAACCCTCTGAGTATGGAACGTT
 TGCTGGTTCCAGGGAAGCCAAGGGACCCCTACTGAATACCCCAAAGCGAGAACGAATGGTGCAT
 GAAGACAGTGGAGGAAAAGAACTTAGAGATTGAGAGGCTTAAGATGAAGCAAAAAGAACTGGAAGCCAAG
 GTGTTGGCCAGGAGGCACCGGATCCAAGGAAAAAGAGAACTCCACAATTCTCCAGCCGCTGCCT
 CTTATAGTGGCACAGTGGCAAAGCCCTGAAAAAGGCTGTAGTATGCCCTGCAACGGATTGAGAAAGCA
 GCGAGAGTCTCAAACCAATCCAGCTCCTGAAGAAAGGCCCAAAGGAAGCTGGAGCCATCCCTGAG
 TCAGAAGCTGTGAAAAGGATGAGGACTACTGGGAGGTACAGATTAGCCAGAGTACTGGCACATGGGC
 GCAAAAAATTATTGGATCTGCTGAATGAAGGCTCAGCTCGGGAAGTGCAGCCTGCAGCGCATTGGCCA
 GAAGAAAGCCAGCTCATCGTGGGCTGGAGGGAGCTCCACGGCCCTTCAAGTGAAGTGGAGGACTTGGAA
 CAGGTGGAAGGCATCTCTGGGAAGCAGGTGGAGTCTTCTGAAGGCGAATCTCTTGAGCTTGGCCGCCA
 GCCAGCACAGCGGTCCCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR209830 representing NM_145588
Red=Cloning site Green=Tags(s)


```
MSLRAKTCPRREMASATSGPGRVCVSKGGLGRRPPLARVRVAVRLRPFMDGETEAKELPCVRAIDSCSLE
VANWKYQETLK YQDFAYGEKSTQQEVVYVGSVQPILRHLLLEGQNASVLAYGPTGAGKTHMLGSPEQPG
VIPRALMDLLQLAREESAEGRPWDVSVAMSYLEIYQEKVLDLLDPASGDLVIREDCRGNILIPGLTQKPI
TSFSDFEQHF LPA SRNRVAVGATRLNQRSSRSHAVLLVKVDQQRERLTPFRQREGKLYLIDLAGESENRRTG
NQGIRLKE SGAIN TSLFVLGKVVDALNQGLPRIPYRDSKLRLLQDSLGGSAHSILIANIAPERRFYQDT
ISALNFTARSKEVINRPF TNE SLQPHALAPVKLSQKELLGPSEAKKAKGPEEESTGSPESTAAPASASQK
LSLLQKLSNMDPAMLENLLSMERLLGSQGSQGTPLLNTPKRERMVLMKTVEEKNEIERLKMKELEAK
VLAQEAPDPREKENTPTILQPPASYS GTVAKPLKKA VVMPLQRIQKQRESSNQIQLLKKGPKRKLEPSPE
SEAVEKDE DYWEVQI SPELLAHGRKLLDLLNEGSARELRSLQRIGQKKAQLIVGWRELHGPFSEVEDLE
QVEGISGKQVESFLKANLLSLAASQHSGPS
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9036_a01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_145588

ORF Size: 1980 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_145588.1](#), [NP_663563.1](#)

RefSeq Size: 2086 bp

RefSeq ORF: 1983 bp

Locus ID: 110033


UniProt ID: [Q3V300](#)

Cytogenetics: 7 69.29 cM

MW: 73.6 kDa

Gene Summary: Kinesin family member that is involved in spindle formation and the movements of chromosomes during mitosis and meiosis. Binds to microtubules and to DNA. Plays a role in congression of laterally attached chromosomes in NDC80-depleted cells.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209830