

## Product datasheet for **MR209830L4V**

### Kif22 (NM\_145588) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Kif22 (NM_145588) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kif22
Synonyms:	AU021460; C81217; Kid; Kif22a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145588
ORF Size:	1980 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209830).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145588.1</a> , <a href="#">NP_663563.1</a>
RefSeq Size:	2086 bp
RefSeq ORF:	1983 bp
Locus ID:	110033
UniProt ID:	<a href="#">Q3V300</a>
Cytogenetics:	7 69.29 cM



[View online »](#)

**Gene Summary:**

Kinesin family member that is involved in spindle formation and the movements of chromosomes during mitosis and meiosis. Binds to microtubules and to DNA. Plays a role in congression of laterally attached chromosomes in NDC80-depleted cells.[UniProtKB/Swiss-Prot Function]