

Product datasheet for **MR209837L3V**

Cxhc1 (NM_028868) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cxhc1 (NM_028868) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cxhc1
Synonyms:	2410002I16Rik; 5830420C16Rik; AI426635; Cfp1; Cgpb; PHF18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028868
ORF Size:	1983 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209837).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028868.1 , NP_083144.1
RefSeq Size:	2537 bp
RefSeq ORF:	1983 bp
Locus ID:	74322
UniProt ID:	Q9CWW7
Cytogenetics:	18 E2
Gene Summary:	Transcriptional activator that exhibits a unique DNA binding specificity for CpG unmethylated motifs with a preference for CpGG.[UniProtKB/Swiss-Prot Function]



[View online »](#)