

## Product datasheet for **MR209840L4V**

### Dbn1 (NM\_019813) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dbn1 (NM_019813) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dbn1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019813
ORF Size:	1983 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209840).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019813.3</a>
RefSeq Size:	2906 bp
RefSeq ORF:	1983 bp
Locus ID:	56320
UniProt ID:	<a href="#">Q9QXS6</a>
Cytogenetics:	13 30.06 cM



[View online »](#)

**Gene Summary:**

Actin cytoskeleton-organizing protein that plays a role in the formation of cell projections (By similarity). Required for actin polymerization at immunological synapses (IS) and for the recruitment of the chemokine receptor CXCR4 to IS (By similarity). Plays a role in dendritic spine morphogenesis and organization, including the localization of the dopamine receptor DRD1 to the dendritic spines (PubMed:25865831). Involved in memory-related synaptic plasticity in the hippocampus (PubMed:25865831).[UniProtKB/Swiss-Prot Function]