

## Product datasheet for **MR209873L4V**

### Sybu (NM\_178765) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sybu (NM_178765) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sybu
Synonyms:	Golsyn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178765
ORF Size:	1992 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209873).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178765.3</a> , <a href="#">NP_848880.2</a>
RefSeq Size:	3152 bp
RefSeq ORF:	1998 bp
Locus ID:	319613
UniProt ID:	<a href="#">Q8BHS8</a>
Cytogenetics:	15 B3.2
Gene Summary:	Part of a kinesin motor-adaptor complex that is critical for the anterograde axonal transport of active zone components and contributes to activity-dependent presynaptic assembly during neuronal development.[UniProtKB/Swiss-Prot Function]



[View online »](#)