

## Product datasheet for **MR209877**

### Clptm1 (NM\_019649) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Clptm1 (NM_019649) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Clptm1
Synonyms:	HS9; N14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR209877 representing NM\_019649  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCGGGCGCAGGAGCGGACGGGGCCGGCAGCGCCGTGGTGGCGGGGGGGAGGCAGCTCTGGTC  
 AGGTGACCAGCAATGGCAGCATTGGGAGGGACAGCCAGCCGAGACCCAGCCTCAGAACCCCTCCACCCCA  
 GCCAGCCCCAACGCCTGGCAGGTCATCAAAGGCGTTCTGTTTCAGGATCTTCATCATCTGGGCCATCAGC  
 AGCTGGTTCGCCGAGGGCCGTCCTCAGGACCAGTCAGGCCCTGGAGGGGGCCCCACGTGTCGCCAGCC  
 GAAACCTATCCCAAAGACTCTAATGAACCTGCACGTGTACATCTCAGAGCAGGAGCACTTCACAGA  
 CTTCAACGCCACCTCGGCGCTCTTCTGGGAGCAGCATGACCTGGTGTATGGTGACTGGACCAGTGGAGAG  
 AACTCTGATGGCTGCTACGAGCACTTTCAGAGCTTGCATCCCACAGAGTGTACAGCAGAATGGCTCCA  
 TCTATATCCATGTCTACTTCAAGAGTGGCTCCATCCAGACCCCGCAGAAGGCCCTGTACCGCCG  
 GCTTGCCACTGTCCACATGTCCCGGATGATCAACAAATACAAGCGTAGGGGATTTAGAAAACCAAGAAC  
 CTAAGTACAGGAGAGACAGAAGCTGACCCAGAAATGATCAAGAGGGCTGAGGACTATGGCCCTGTAGAAG  
 TAATCTCGCACTGGCACCCCAACATCACCATCAACATCGTGGATGACCACACACCATGGGTGAAGGGCAG  
 TGTGCCACCCCACTGGACCAATATGTGAAGTTTGTGCGGTGAGTGGTGACTACTACCCCATCATCTAT  
 TTCAATGACTACTGGAACCTGCAGAAGGATTACTACCCATTAATGAGAGTCTGGCCAGCCTACCACTGC  
 GTGTCTCTTTCTGCCGCTGTCACTATGGCGCTGGCAACTCTATGCTGCACAGAGTACCAAGTCAACCTG  
 GAATCTTAGGTGACGAGCTGTATGAGCAATCGGATGAGGAGCAGGATTCAGTGAAGTTGCTCTGCTG  
 GAAACCCAGCCGTATCTTCTGGCACTCACCATCATCGTCTCCATTGTCCACAGTGTCTTTGAGTTCCTGG  
 CTTTCAAAAATGACATTCAGTTCTGGAACAGTCGTCAGTCTTGGAGGGGCTGTCAGTGGCTCCGCTCTT  
 CTTTGGCGTTTTCCAGTCCTTTGTGGTTCTCCTTTACATCCTGGACAATGAGACCAACTTTGTGGTACAG  
 GTCAGCGTCTTCATTGGGGTCTCATCGACCTCTGGAAGATTACCAAAGTCATGGACGTCGGCTGGACC  
 GGGAGCATAGAGTGGCAGGAATCTTCCCTTGCCCAACCTCAAAGACAAGTCCACCTATATTGAATCTTC  
 AACCAAAGTATATGATGATATGGCGTTCGGTACCTGTGATGGATCCTCTTTCCCTTGGTGGGCTGCTAT  
 GCCGTCTACAGCCTTCTGTACCTGGAGCACAAGGGCTGGTACTCCTGGGTGCTCAGCATGCTCTACGGTT  
 TCCTGCTGACCTTTGGCTTCATCACTATGACACCACAGCTCTTCATCACTATAAGCTCAAGTCTGTGGC  
 CCATTTGCCCTGGCGGATGCTCACCTACAAGGCCCTCAACACTTTCATTGATGACTTATTTGCCTTTGTC  
 ATCAAGATGCCTGTGATGATACCGGATCGGCTGCCTGCGGGATGATGTGGTTTTTTCATCTACCTCTACC  
 AACGGTGGATCTACCGCGTCGACCCACCCGTGTGAACGAGTTCCGGTATGAGTGGAGAGGATGTCTCAGC  
 AGCGGCATCCAGGGCCAGGCATCCACTGCAGCAGGGGCTCTCACACCTGCACCCAGCACAGCTGTGAGC  
 GGAGAGGATGCCTCCACAGTCCCAAAGGCCACTTCTGGGGCTGCACTGCCAGCCAGCCCCAGGAAGCCC  
 CTCCAAAGCCAGCAGAGGACAAGAAAAGGAT

**ACGCGT**ACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR209877 representing NM\_019649  
Red=Cloning site Green=Tags(s)


MAAAQEADGAGSAVVAAGGGSSGQVTSNGSIGRDTPAETQPQNPQPAPNAWQVIKGVLFRIIFIWAIS  
 SWFRRGSPQDQSGPGGAPRVASRNLFPKDTLMNLHVYI SEHEHFTDFNATSALFWEQHDLVYGDWTSGE  
 NSDGCYEHFAELDIPQSVQNGSIYIHVYFTKSGFHPDPRQKAL YRRLATVHMSRMINKYKRRRFQKTKN  
 LLTGETEADPEMIKRAEDYGPVEVISHWHPNITINIVDDHTPWVKGSVPPPLDQYVKFDAVSGDYPIIY  
 FNDYWNLQKDYYPINESLASLPLRVSFCLSLWRWQLYAAQSTKSPWNFLGDEL YEQSDEEQDSVKVALL  
 ETSPYLLALTIIVSVHVSVEFLAFKNDIQFWNSRQSLEGLSVRSVFFGVFQSFVLLYILDNETNFVVQ  
 VSVFIGVLIDLWKITKVM DVRLDREHRVAGIFPCPTFKDKSTYIESSTKVYDDMAFRYLSWILFPLLGCY  
 AVYLLYLEHKGWYSWLSMLYGFLLTFGFITMTPQLFINYKLSVAHL PWRMLTYKALNTFIDDLFAFV  
 IKMPVMYRIGCLRDDVFFIYLQYRWIYRVDPTRVNEFGMSGEDVSAASRAQASTAAGALTPAPSTAVS  
 GEDASTVPKATSGACTASQPQEAPPKPAEDKKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9014\\_a04.zip](https://cdn.origene.com/chromatograms/mm9014_a04.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_019649

**ORF Size:** 1992 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_019649.2](#), [NP\\_062623.2](#)

**RefSeq Size:** 4147 bp

**RefSeq ORF:** 1995 bp

**Locus ID:** 56457


**UniProt ID:** [Q8VBZ3](#)

**Cytogenetics:** 7 A3

**MW:** 75.7 kDa

**Gene Summary:** May play a role in T-cell development.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR209877