

## Product datasheet for **MR209887**

### Zfp760 (NM\_001008501) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp760 (NM_001008501) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zfp760
Synonyms:	9430018I06; BC029103
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR209887 representing NM\_001008501  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGCTGCTTCTCCAGGTCTGTTGACATTCGGGGATGTGGCTTTGGACTTCTCCAGGAGGAGTGGGAAT  
GCCTGGACTCTGCTCAGAGGGCTCTGTACATTGATGTGATGTTGGAAAATTACATCAACTTGGTCTCTGT  
TGAGAAGCATCACATGTGTGGTAAATATGAGAAGTCTTGGATCAACACTCAGAGCATCAACATGTGAAT  
ATCCAAGAAAAGTCTTACAAATATAAAGAGCTTGGCAAAATGGTTCATGAATCCTCTCCATGTACACCTT  
ATAACACAAGTAGTACTGCAGAAAAGTGAACAAGTACAGCTACCAGGATGCTTCCATTGAATCATCCCT  
AATCAGACATAATAGTGTACACACTGGAGAAGAACTTTTCAAATATAAAGAGTGTGGCAAAATGTTTATGT  
TCCAACATTTGTCCAATCAGAGTATCCACGATATAAAGATGGAACACAAGTATACAGAAAGTGATAAGC  
TTTGTGGCTCTACACACAGACTTGTACAGCAAAGAACTCATAGTGGAGAGAAAGCACACAATTGTGGAA  
ATGCTTTAGGAGTTGCTCAAGCCTCAGTACACATCAGAGACTTCACACTGGAGAGAAACCATGGAATGT  
AAAGAATGTGGTAAGTCCCTTAAAGTGGTTGTCTCATCTTAAAGCCATTGCAGAAATCCACTGGAGAAA  
AACCTTTCAAATGTAATGAATGTGATAAAGGCTTTTACGTTATCTAGACTTAAAAGTCATTATAGAAC  
CCACTACTGGAGAGAAGCCTTACAAATGTGAGGAATGCGGTAAGTCTTTTCTCGTTGTCAAAGCTTAAA  
AACCATTATAGAATCCATACAGGAGAGAAAACCTTACAAGTGAAGGAATGTGGCAAAATCTTTTACCATT  
GCACGGTCTTAGTACACATCAGAAAATTCATACTGGGGAGAAAACCTCACAAATGCAAAAAGTGGTAA  
ATCCTTTTCTGCTGTCAACTCTTAAAGAACATTTTCAAGTCCACACTAGAGAGAAGCCTTACAAGTGC  
AAAGAATGTGGTAAGTCTTCCATTGGTTGTGAGGCTTAAAAGCCACTACAGAATCCACTGGAGAGA  
CACCTTACAAATGTAAGGAGTGTGACAAGTCAATTTCTCAGTTATCAAGACTTAAAAGCCATTATAGAAC  
CCACTACTGGAGAGAAAATCTTCAAATGTTTAAAGTGTGACAATACTTTACCCAATGCACAGCTCTTAAA  
AGACATCAGAAAATTCATACTGATGAAAAACCTACAAATGTGAAGAGTGCAGCAAAATGTTTTTATGTGT  
TGTCTCACCTAAAAGCCATTGCAGAACCATACTGGAGAGAAGCCTTACAAATGTGAGGAATGTGGTAA  
GTGCCTTTCTCAGTTGTCAAACCTTAAAAGCCATTATAGAATCCATACAGGAGAGAAAACCTTATAAGTGT  
AAGGAATGCGACAAAATCCTTTACCATTGCACAGCTCTTAGTACACATCAGAAAATTCATACTGGAGAGA  
AACCTACTCATGCAAAGACTGTGGTAAATCCTTTTATCTGTTATCAAGACTTAAAGAACATTACAGAGT  
CCACACTGGAGAGAAAACCGTACAATTGTAAGGAATGTGGTAAGTCTTCCATTGGTTGTCAAACCTTAAA  
AGACATTCAGAAATTCACACTGGAGAGAAAACCTTACAAATGTGAGGAATGTGACAATCCTTTACCCT  
GCTCAAGTCTTAAAGCTCATCAGAAAATTCATGCTGGGGTAAAACCTACAAATGTGAAGAATGTAGCAA  
GTCCTTTTATTGGTTGTCAAACCTTAAAACCCATTCCAGAATCCACACTGGAGAAAAACCTTACAAATGT  
GGTAAATGTGACAATCCTTTACCAGCGCTCTCTTCTTACAAAGCACACAGGATCCACTACTAGAGAAA  
AACATGAAAACAGGAAATATGTGAGAATTTTGATGAA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR209887 representing NM\_001008501  
Red=Cloning site Green=Tags(s)


```

MAASPLLLTFGDVALDFSQEEWECLDSAQRALYIDVMLENYINLVSVEKHHMCGKYEKVLDQHSEHQHVN
IQEKSYKYKELGKMVHESPCTPYNTSSTAENCNKYSYQDASIESSLIRHNSVHTGEELFKYKECGKCLC
SNICPNQSIHDIKMEHKYTESDKLCGSTHRLVQQRTHSGEKAHNCCKFRSCSSLSTHQRLHTGEKPWKC
KECGKSFKWLHLKAHCRIHTGEKPFKCNEDKGF SRLSRLKSHYRTHTGEKPYKCEECGKCF SRLSRLK
NHYRIHTGEKPYKCEGKSFTHCTGLSTHQIHTGEKPHKCKNCGKSFHLLSTLKEHFRVHTREKPYKC
KECGKFFHWL SGLKSHYRIHTGETPYKCKE CDKSFSQLSRLKSHYRTHTGEKSCKLECDKYFTQCTALK
RHQKIHTDEKPYKCEEC SKFYVL SHLKSHCRTHTGEKPYKCEECGKCLSQLSKL SHYRIHTGEKPYKC
KECDKSFTHCTALSTHQIHTGEKPYKCKDCGKSFYLLSRLKEHYRVHTGEKPYNCKE GKFHWL SNLK
RHSRIHTGEKPYKCECDKSFTHCSSLKAHQKI HAGVKPYKCEEC SKSFYWL SNLKTHSRIHTGEKPYKC
GKCDKSF TQRSLLT KHNRIHTREKHENRKYVRTFDE
    
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001008501

**ORF Size:** 1998 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001008501.2](#), [NP\\_001008501.2](#)

**RefSeq Size:** 4209 bp


**RefSeq ORF:** 2001 bp

**Locus ID:** 240034

**Cytogenetics:** 17 A3.2

**MW:** 77.9 kDa

**Product images:**



Circular map for MR209887