

Product datasheet for **MR209914L3V**

Dym (NM_027727) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dym (NM_027727) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dym
Synonyms:	1810041M12Rik; 4933427L07Rik; C030019K18Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027727
ORF Size:	2010 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209914).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027727.2
RefSeq Size:	2471 bp
RefSeq ORF:	2010 bp
Locus ID:	69190
UniProt ID:	Q8CHY3
Cytogenetics:	18 E2-E3
Gene Summary:	Necessary for correct organization of Golgi apparatus. Involved in bone development. [UniProtKB/Swiss-Prot Function]



[View online »](#)