

Product datasheet for **MR209922L3V**

Zbtb5 (BC057545) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zbtb5 (BC057545) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zbtb5
Synonyms:	5930421110; 9430083K24Rik; AI646847; mKIAA0354
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC057545
ORF Size:	2010 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209922).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC057545 , AAH57545
RefSeq Size:	3871 bp
RefSeq ORF:	2012 bp
Locus ID:	230119
Cytogenetics:	4 B1
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]


[View online »](#)