

Product datasheet for **MR209925L3V**

Tbc1d15 (NM_025706) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tbc1d15 (NM_025706) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tbc1d15
Synonyms:	4432405K22Rik; AU022358; AU023730; AW557774; Ly6dl; Rab7-GAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025706
ORF Size:	2013 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209925).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025706.3 , NP_079982.3
RefSeq Size:	3779 bp
RefSeq ORF:	2016 bp
Locus ID:	66687
UniProt ID:	Q9CXF4
Cytogenetics:	10 D2
Gene Summary:	Acts as a GTPase activating protein for RAB7A. Does not act on RAB4, RAB5 or RAB6. [UniProtKB/Swiss-Prot Function]



[View online »](#)