

Product datasheet for **MR209926**

Spg20 (NM_144895) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spg20 (NM_144895) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Spg20
Synonyms:	A1840044; C79168; mKIAA0610; SPARTIN; TAHCCP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR209926 representing NM_144895
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGAGAGAGAACCAGAAAATGGAGAACCTGCTGAAATTAAGATCATCAAGGAAGCATACGAGAAGGCCT
 TTATGTTTGTAAATAAAGGACTGAACACGGACGAGTTAGGGCAGAAGGAAGAAGCAAAGAACTATTACAA
 GCAAGGAATAGGACACCTGCTTAGGGGCATCAGCATCGCGCTGCAGAGCCGGGCACACAGGCCCTGCG
 TGGGAGGCTGCCAGGCAGATGCAACAGAAAATGAAGGAAACGCTTCAGAACGTCCGCACCCAGGCTAGAGA
 TCCTGGAGAAGGGCCTTGCCACTTCTACGGAATGACCTTCAGGATGTGCCAAGCTATACCCAGAGTT
 CCCGCCTAAAGACGCGTGAAGAAGTCCCAGAGCAGGAGTCCGTTAGTACAGCGCCCCAGCGTGTGAA
 GTCGATGGGAGCGCTCAGCCGCGTGTGCCGGCCAAGTGGTGCCTTCTGCTCTGCCGTACCATCGC
 CCAGCTGCCCTGCCGAAGCGCCACCAGCCTACAGCCACAGGCTGCTGAGGGCCACTACACTGTTTCCTA
 TGGAACAGATTCCGGGGAGTTTTTCATCCGTGGGAGAAGACTTCTACAGGAACCGTCCCAGCCGCCCT
 CTTGAGACCTTGGGGCTCGATGCAGATGAGTTGATCTTGATACCAAACGGAGTGCAGATCTTTTTGTAA
 ACCCTGCGGGAGAGGTCAGTGCTCCCTCGTACCCTGGGTATCTGCGGATCGTGCCTTCTTGATAATTC
 TCTGGATACAGTTCTGAACCGCCTCCTGGGTTTCTCCAGGTGTGTGACTGGTTATACCTCTAGTTCT
 GATAGTCAACAGTTCTGAAGTGCAGTGTGGGAGCCTACATGTTCCAGATACCATGTTACAAGCTGCTG
 GGTGCTTCGTAGGGTTGTTCTTTCATCTGAGCTCCCAGAAGTACAGAGAGCTGTTTGGAGTCTGTT
 GAGACAGATGTCTGACCTTCGGCTCCAGGCCAACTGGAATAGAGAAGAAGATGAATCCAAATCCCTGGG
 AGATCAAGTCAACCTCTGAACCACCGAAGGATGCCCGCCATAAAGGAAAGCGTGGAAAAAGACCAAAGATTC
 GTTCTCAATAGACCAAGGCAGCAAGGATGCCCGCCATAAAGGAAAGCGTGGAAAAAGACCAAAGATTC
 AAGTGAAGAAGTTAATCTGAGCCAGATTGTGCCCTGTGAGCCGAGTTCAGAAGAAAAATCGAAAGAATTG
 CCTGAATGGAGCGAGAAAGTGGCCACAACTCCTGTCAGGTGCTTCTGGGTGAGCTGGGGTTTGTCA
 AAGGTGCTGAATTTACTGGCAAAGCAATTGAGAAAGGTGCTTCTAAACTCAGAGAGCGCATTGAGCCGGA
 GGAAAAGCCAGTGGAGGTGAGCCAGCTGTACCAGGGGTCTCTACATAGCGAAGCAGGCAACCGGAGGA
 GCAGCCAAAGTCAAGCAGCTCCTGGTTGACGGCTCTGCACAGTAGCGAATTGTGTTGGTAAAGAACTG
 CTCCACACGTCAAGAAGCACGGAAGCAAACCTGGTTCCGGAATCCCTGAAGAGGGACAAAGACGGGAAATC
 TGCCCTGGATGGCGCCATGGTCGTGGCAGCGAGTAGTGTCAAGGATTCTCAACTGTCTGGCAAGGGTTG
 GAATGTGCAGCTAAATGCATCGTTAACAACGTGTCTGCAGAACTGTACAGACTGTGAGATACAAGTATG
 GGCACAACGCAGGAGAAGCTACACATAACGCAGTGGATTGAGCTATCAATGTTGGGCTCACTGCCTACAA
 TATTGACAACATCGGCATCAAAGCCATGGTGAAGAAAAGTCCAAAGCAGACGGGACACACTCTCCTCGAG
 GACTACCAGATCGTTGAGCGCCCTCAGAGGGAGAGTCAAGGAGGAGCCACAAGCAGAGAGGAGAAAGAG
 ACATAGGAAAGCAGGTGGAGGAAGAAAGCCAGGGGCAGGAAAGAAAGGACAAA

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR209926 representing NM_144895
 Red=Cloning site Green=Tags(s)


MEREPENGEPAEIKIIEAYEKAFMFVNKGLNTDELGQKKEAKNYYKQGIGHLLRGISIAAAEPGHTGPA
 WEAARQMQQMKETLQNVTRLEILEKGLATSLRNDLQDVPKLYPEFPPKDACCKSPEQESVSTAPQRAE
 VDGASAAACAGPSGAPSALPVPSPSCPAEAPPAYSPQAAEGHYTVSYGTDSEFSSVGEDFYRNRSQPPP
 LETLGLDADELILIPNGVQIFFVNPAGEVSAPSYPGYLRIVRFLDNLSDTVLNRPPGFLQVCDWLYPLVP
 DRSPVLKCTVGAYMFPDMLQAAGCFVGVVLSSELPEDDRELFEDLLRQMSDLRLQANWNREDEFQIPG
 RSSHPSEPPKEASGTDVRQSSSSGSSIDQGSKDARHKGKRGKTKDSSEEVNLSQIVPCEPSSEEKSKEL
 PEWSEKVAHNILSGASWVSWGLVKGAEFTGKAIQKGASKLRERIQPEEKPVEVSPAVTRGLYIAKQATGG
 AAKVSQLLVGCTVANCVGKELAPHVKKHGSKLVPESLKRDKDGKSALDGAMVVAASSVQGFSTVWQGL
 ECAAKCIVNNVSAETVQTVRYKYGHNAGEATHNAVDSAINVGLTAYNIDNIGIKAMVKKAKQTGHTLLE
 DYQIVERPQRESQGGATSTEGRRDIGKQVEEEKPGAGKKDK

SGP TRTRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-RsrII

Cloning Scheme:



ACCN: NM_144895

ORF Size: 2013 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_144895.2](#), [NP_659144.1](#)

RefSeq Size: 3797 bp

RefSeq ORF: 2016 bp

Locus ID: 229285


UniProt ID: [Q8R1X6](#)

Cytogenetics: 3 C

MW: 73.1 kDa

Gene Summary: May be implicated in endosomal trafficking, or microtubule dynamics, or both. Participates in cytokinesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209926