

Product datasheet for **MR209932L3V**

Aacs (NM_030210) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aacs (NM_030210) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aacs
Synonyms:	2210408B16Rik; SUR5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030210
ORF Size:	2016 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209932).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030210.1
RefSeq Size:	3175 bp
RefSeq ORF:	2019 bp
Locus ID:	78894
UniProt ID:	Q9D2R0
Cytogenetics:	5 G1.1
Gene Summary:	Activates acetoacetate to acetoacetyl-CoA. May be involved in utilizing ketone body for the fatty acid-synthesis during adipose tissue development (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)