

Product datasheet for **MR209940**

Raver2 (NM_183024) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Raver2 (NM_183024) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Raver2
Synonyms:	A430091O22Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR209940 representing NM_183024
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGGCGGCGGGGAGGCGGCGGGGGCGGGGCTCGGATCCGGACCCTCTGCTGGGACGGCGGGG
 AGGCCGCGGAGCCCGCGCTGCGCCCGGGGAGGTGGCTGCGCTGCACCCAGGAGGTGGCTGCACGGCT
 GCAGCGGATGCGACGGGAGCTGAGCAACCGCGGAAAAATCCTAGTGAAAAACCTGCCCCAGGACAGCAGC
 TCCAGGAGGTTTCATGAGTTATTACAGGACTATGAGCTAAAGTATTGCTATGTGGACAGGAATAAACGAA
 CAGCTTTTGTACCTTATTGAATGGGAGCAGGCGCAGAGCGCGATTGAGAGTTTCATCAGTTTCTTT
 CCGAGGCAGAGAGTTGACAGTGCAGCTGCAGCCACAGACGCTCTCTGTGCATACCAACCTGCCATC
 TCCTTCACCTTAGAAGAGTTTGAAGAACTCGTCCGTGCCTATGGCAACATCGAGAGATGTTTTTGGTCT
 ACAGTGAAGTTACTGGCCACTCAAAGGCTACGGGTTTGTGGAGTACATGAAGAAGGACTTTGCTGCCAA
 GGCCAGACTGGAGCTACTGGGCGGAGATGGGAGCCTCGGCGCTCTTCGCACAGTGGATGGATGTTAAT
 CTCTTGGCCTCGAGCTCATCCATTCTAAGTGCCTTTGTATAGATAAACTCCCCAGTGACTACAGCGATT
 CGGAAGAGCTCTTGCAACTTTTTCCGGTATCCATAAACCTGTGTTTTGCCAGCTTGACAGGATGAAGG
 AAGTCACGGTGGTGGCTTTGCGGTGGTTGAATATAGCACTGCGGAGCATGCCAAGAGGTTGAGCAGGTG
 GCAGATGGCATCACCATCAAGGGGAGCCAAAGTTCAACTGTCTTCTGCGCCCGAGGACACAGGGCGAA
 GCACTCTAGCAGTCTCATAGCAGCTCAGCGCGGATGCACAGTAATCAAAGGGATTACTGCCAGAACC
 CAATCCGGTACAGATTATGAAAAGTTTAAACAACCTGCCATGTTACAAGTTCTTCTGCAGCCCCAGCTA
 TCGGACGAGCCATGAAGCCAGTTCTCGGAGTGCAGCCAGCTTATCTCACCTGCTGAGTCCGTCAGTCC
 CCTCCGCCATCTTACACTTCAGTAAAGCACAGCAGAGCTCAGCTGTGGTAATACATCGTCTTAATCCT
 TCAGAACCTGTCTCCCTTGCCCTTGATACAGCAACAGTTGATGAAGTTTGACAATGCCACACTAACAAAT
 AAGCCTGGCTTGCTGGGAGAGCCTCCTGCCATGGTACTGCAGCCTGCTCTGGCGATAGGGCCACCGCTTC
 CACTGAAAACGGACCTGGGACACCATGGCGAAGCACACAAAACATCTAACCTGATTCCACCTCAGACGAC
 ATTGGCAGCCGGGCTGGGCATGCTGCCGTTCTTCTCAAACCAGCTCCCTGCTGGACAAGCCGGCCGGGG
 CGTGGGACCACTCAGGAGAAACAGTCCGCTTCGGTGAGCATTAGCGAAGCTAGCTTCTCAGGGTCGCAGC
 ATTACCTGCAGACCTTCCAGGCCTTCTGCTGGAGGCCCGCTGACGGGCAACCAGAAGACACCGCAGAG
 CCAGCCAAAAGGCACAGAGGTTGCCTCTAAGAATCAGACTTCACTCTTGGGAGAACCACAAAAGAGATC
 CGGCTCAGTAAAAATCCATATTTGAATTTGGCAAGTGTGCTACCCAGTGTGTGCCTATCCACTGCAGGTA
 AAGGCATGCCTCCGAAGACTGGAATTGCGAGCAACATTCTGGATGCGATCTCTCAGGGAAGTGAGTCCCA
 GCATGCACTGGAGAAGTGCATCGCTATTCTCCATCCATCGAGGATTATGCCAGGCATCTCTTGGAGG
 AATGAGAAGCGAGGCTCCTCTATCTCTGCCCCAGAAGGAGGCCCGGTAGAACTCGCTGGCCAGC
 ACCCTCAGGACACCGGAGTGAGCTACACAGAACTACTTGAAAAAGAAGCGTGTGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR209940 representing NM_183024
Red=Cloning site Green=Tags(s)


MAARGGGAGGAGSGSGPSAGTAGEAAEPALRPGEVAALHPQEVAARLRMRRELSNRKILVKNLPQDSS
SQEVHELLQDYELKYCYVDRNKRTAFVTLNNGEQASAIQRHFQFSFRGRELTVQLQPTDALLCITNLP
SFTLEEFELVRAYGNIERCFLVYSEVTGHSKGYGFVEYMKKDFAAKARLELLGRQMGSALFAQWMDVN
LLASELIHSKCLCIDKLPSDYSDEELLQLFSGIHPVFCQLAQDESGHGGFVAVEYSTAEHAEVQVQV
ADGITIKGSQVQLSFCAPGAPGRSTLAVLIAAQRAMHSNQKGLLPEPNPVQIMKSLNPNAMQLVLLQPQL
CGRAMKPVLGVPAPSLSHLLSPSLPSAILHFSKAQQSSAVGNTSSLILQNLSPPLLIQQQLMKFDNAHTNN
KPGLLGEPPAMVLQPALAIGPPLPKTDLGHHGEAHKTSNLIIPPQTTLAAGLGMLPFFSNQLPAGQAGPG
RGTTQEKKQSASVSISEASFSGSQHYLQTFPGLPAGGPLTGNQKTPQSQPKGTEVASKNQTSLLEGPKKI
RLSKNPYLNLASVLPVCLSTAGKGMPKPTGIASNILDAISQGSSEQHALEKCIAYSPSIEDYAQASSLR
NEKRGSYLSIAPEGGPVLAGOHPDGTGVSYTETYLKKKRVY

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9104_c07.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_183024

ORF Size: 2019 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_183024.1](#), [NP_898845.1](#)

RefSeq Size: 4006 bp

RefSeq ORF: 2022 bp

Locus ID: 242570


UniProt ID: [Q7TPD6](#)

Cytogenetics: 4 C6

MW: 72 kDa

Gene Summary: May bind single-stranded nucleic acids.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209940