

Product datasheet for **MR209965L3V**

Nub1 (NM_016736) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nub1 (NM_016736) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nub1
Synonyms:	4931404D21Rik; 6330412F12Rik; BS4; NY-REN-18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016736
ORF Size:	1845 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209965).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016736.2 , NP_058016.2
RefSeq Size:	3455 bp
RefSeq ORF:	1845 bp
Locus ID:	53312
UniProt ID:	P54729
Cytogenetics:	5 A3
Gene Summary:	Specific down-regulator of the NEDD8 conjugation system. Recruits NEDD8, UBD, and their conjugates to the proteasome for degradation (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)