

Product datasheet for **MR210002L3V**

Gfpt1 (NM_013528) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gfpt1 (NM_013528) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gfpt1
Synonyms:	2810423A18Rik; AI324119; AI449986; GFA; GFAT; GFAT1; GFAT1m; Gfpt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013528
ORF Size:	2046 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210002).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013528.2 , NP_038556.1
RefSeq Size:	6248 bp
RefSeq ORF:	2046 bp
Locus ID:	14583
UniProt ID:	P47856
Cytogenetics:	6 37.81 cM
Gene Summary:	Controls the flux of glucose into the hexosamine pathway. Most likely involved in regulating the availability of precursors for N- and O-linked glycosylation of proteins. Regulates the circadian expression of clock genes ARNTL/BMAL1 and CRY1.[UniProtKB/Swiss-Prot Function]



[View online »](#)