

## Product datasheet for MR210008L3

### Dtna (NM\_207650) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Dtna (NM_207650) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Dtna
Synonyms:	2210407P21Rik; a-DB-1; adbn; Dtn; DTN-A; Gm19389
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210008).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_207650
ORF Size:	2049 bp



[View online »](#)

