

Product datasheet for **MR210021L4V**

Slco2b1 (NM_175316) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slco2b1 (NM_175316) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slco2b1
Synonyms:	AI060904; AI852653; OATP-B; Slc21a9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175316
ORF Size:	2049 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_175316.3 , NP_780525.1
RefSeq Size:	4381 bp
RefSeq ORF:	2052 bp
Locus ID:	101488
UniProt ID:	Q8BXB6
Cytogenetics:	7 E1
Gene Summary:	Mediates the Na(+)-independent transport of organic anions such as taurocholate, the prostaglandins PGD2, PGE1, PGE2, leukotriene C4, thromboxane B2 and iloprost. [UniProtKB/Swiss-Prot Function]



[View online »](#)