

Product datasheet for **MR210037L3V**

Zfp12 (BC053080) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zfp12 (BC053080) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zfp12
Synonyms:	C530015C18; Krox-7; Zfp-12; Znf12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC053080
ORF Size:	2058 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210037).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC053080.1
RefSeq Size:	4315 bp
RefSeq ORF:	2060 bp
Locus ID:	231866
Cytogenetics:	5 G2
Gene Summary:	Transcriptional repressor which suppresses activation protein 1 (AP-1)- and serum response element (SRE)-mediated transcriptional activity.[UniProtKB/Swiss-Prot Function]


[View online »](#)