

Product datasheet for **MR210060L3V**

Pgghg (NM_145387) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pgghg (NM_145387) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pgghg
Synonyms:	5730511L01Rik; Ath11; BC023151
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145387
ORF Size:	2073 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210060).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145387.2
RefSeq Size:	2975 bp
RefSeq ORF:	2073 bp
Locus ID:	212974
UniProt ID:	Q8BP56
Cytogenetics:	7 F5
Gene Summary:	Catalyzes the hydrolysis of glucose from the disaccharide unit linked to hydroxylysine residues of collagen and collagen-like proteins.[UniProtKB/Swiss-Prot Function]



[View online »](#)