

Product datasheet for **MR210085L3V**

C1s (BC018319) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C1s (BC018319) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	C1s
Synonyms:	MGC37773, C1sa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC018319
ORF Size:	2082 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210085).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC018319.1
RefSeq Size:	2855 bp
RefSeq ORF:	2084 bp
Locus ID:	50908
Cytogenetics:	6 F2
Gene Summary:	C1s B chain is a serine protease that combines with C1q and C1r to form C1, the first component of the classical pathway of the complement system. C1r activates C1s so that it can, in turn, activate C2 and C4 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)