

Product datasheet for MR210117L3

Acs11 (NM_007981) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Acs11 (NM_007981) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Acs11
Synonyms:	Ac; Acas; Acas1; Acs; Fac; Fac12; FACS; LACS 1; LACS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210117).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_007981
ORF Size:	2097 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_007981.3 , NP_032007.2
RefSeq Size:	3891 bp
RefSeq ORF:	2100 bp
Locus ID:	14081
UniProt ID:	P41216
Cytogenetics:	8 B1.1
Gene Summary:	The protein encoded by this gene belongs to a family of acyl coenzyme A synthetase proteins, which convert long chain fatty acids to acyl CoA products via an ATP-dependent pathway. This enzyme is enriched in heart, liver and adipose tissue, where it functions in lipid synthesis and mitochondrial and peroxisomal beta-oxidation. In addition, it is expressed in monocytes and macrophages where it appears to have a functionally distinct role in mediating inflammatory and innate immune responses. A pseudogene of this gene is found on chromosome 5. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]

Product images:



Circular map for MR210117L3