

## Product datasheet for **MR210119L3V**

### Adam5 (BC094237) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adam5 (BC094237) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adam5
Synonyms:	tMDC II; tMDCII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC094237
ORF Size:	2097 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210119).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC094237</a> , <a href="#">AAH94237</a>
RefSeq Size:	2290 bp
RefSeq ORF:	2099 bp
Locus ID:	11499
Cytogenetics:	8 13.26 cM



[View online »](#)

**Gene Summary:**

This gene encodes a member of a disintegrin and metalloprotease (ADAM) family of endoproteases that play important roles in various biological processes including cell signaling, adhesion and migration. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional protein. This gene is located in a cluster of related ADAM genes on chromosome 8. Alternative splicing results in multiple transcript variants encoding different isoforms, some of which may undergo similar processing. [provided by RefSeq, May 2016]