

## Product datasheet for **MR210122L4V**

### **Dnm1l (NM\_001025947) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Dnm1l (NM_001025947) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dnm1l
Synonyms:	6330417M19Rik; AI450666; Dlp1; Dnm; Dnmlp1; Dr; Drp1; pyt; python
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001025947
ORF Size:	2100 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210122).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001025947.1</a> , <a href="#">NP_001021118.1</a>
RefSeq Size:	4034 bp
RefSeq ORF:	2100 bp
Locus ID:	74006
UniProt ID:	<a href="#">Q8K1M6</a>
Cytogenetics:	16 A2



[View online »](#)

**Gene Summary:**

This gene encodes a member of the dynamin family. The encoded protein is localized to the cytoplasm and mitochondrial membrane, is involved in mitochondrial and peroxisomal division, and is essential for mitochondrial fission. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 2. [provided by RefSeq, Feb 2013]