

Product datasheet for **MR210127L4V**

Capn2 (NM_009794) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Capn2 (NM_009794) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Capn2
Synonyms:	AI326419; Capa-2; Capa2; m-calpain; m-calpin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009794
ORF Size:	2100 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210127).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009794.3
RefSeq Size:	3205 bp
RefSeq ORF:	2103 bp
Locus ID:	12334
UniProt ID:	O08529
Cytogenetics:	1 H5



[View online »](#)

Gene Summary:

Calcium-regulated non-lysosomal thiol-protease which catalyzes limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction. Proteolytically cleaves MYOC at 'Arg-226' (By similarity). Proteolytically cleaves CPEB3 following neuronal stimulation which abolishes CPEB3 translational repressor activity, leading to translation of CPEB3 target mRNAs (PubMed:22711986).[UniProtKB/Swiss-Prot Function]