

Product datasheet for **MR210153L4V**

Gys2 (NM_145572) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gys2 (NM_145572) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gys2
Synonyms:	BC021322; LGS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145572
ORF Size:	2112 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210153).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145572.2 , NP_663547.2
RefSeq Size:	2671 bp
RefSeq ORF:	2115 bp
Locus ID:	232493
UniProt ID:	Q8VCB3
Cytogenetics:	6 G2
Gene Summary:	Transfers the glycosyl residue from UDP-Glc to the non-reducing end of alpha-1,4-glucan. [UniProtKB/Swiss-Prot Function]



[View online »](#)