


## Product datasheet for MR210176L4

### C1r (BC004637) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	C1r (BC004637) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	C1r
Synonyms:	C1rb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210176).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	BC004637
ORF Size:	2121 bp



[View online >](#)

