

Product datasheet for **MR210179**

Eomes (NM_010136) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eomes (NM_010136) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Eomes
Synonyms:	C77258; TBR-2; Tbr2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR210179 representing NM_010136
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCACTTGGGAGAGCAGCTCCTGGTGAGCTCGGTGAACCTGCCCGCGCGCACTTCTACTCGTGGAGA
 GTGCTCGCGGAGGAGGAGGAGGAGGAGGCGGAGGAGGAGGGGAGGCGGAGGGAGCGTCAGCCTCCTCCC
 CGGTGCTGCCCCCTCGCCCCAGAGGCTGGACTTAGACAAAGCGTCCAAGAAGTTTCCGGGAGTCTCCCG
 TGCCAGGCGGGGAGCGCAGAACCCGAGGCGCCGGCGCGGGGCCCCCGCGGCCATGCTCAGTGACGCGG
 ACGCTGGGACACCTTCGGGAGCACCTCGGCGGTGGCCAAGCCCGGGCCCCCGGACGCGCCGAAGGGCTC
 CCCGTGCGCGGAGGAGGAGCTGCCCTCCGCGCCACCGCCGCGGCCACCGCGCGTACTCCATGGACAGC
 CTGAGCTCCGAGCGCTACTACCTCCCGTCGCGGGACCGCAGGGCTCCGAGCTCGCCGCGCCCTGCTCGC
 TCTTCCAGTACCCGGCGCGGCCGGAGCAGCCACGGACCCGTGTACCCCGCTCCAATGGCGCGCGCTA
 CCCCTACGGCTCCATGCTGCCCCCGGTGGATTCCCCCGCGCGGTGTGCCCGCCCGAGGGCGCAGTTC
 GGCCCCGCTGCGGGTTCGGGGAGCGGCGTGGTAGCAGCGCGGTGGTGC CGGCGGTCTGGCGCCTATC
 CCTACGGCCAGGGTTCTCCGCTCTACGGGCCATACGCCGGAACCTCAGCGGCCGGGTCTTGTGGAGGATT
 GGGGGGCTAGGGGTGCCTGGCTCCGGCTTCCGCGCCACGTCTACCTGTGCAACCGGCCCTATGGCTC
 AAATCCACCGGCACAAACTGAGATGATCATCAACAAACAGGGCAGGCGCATGTTTCTTTCTTGGAGCT
 TCAACATAAACGGACTCAACCCACCGCCACTACAATGTTTTCGTGGAAGTGGTCTGGCCGACCTAA
 CCACTGGCGCTTCCAGGGGGCAAGTGGGTGACCTGCGGCAAAGCGGACAATAACATGCAGGGCAATAAG
 ATGTACGTTACCCAGAATCTCCTAACACTGGCTCCCCTGGATGAGGCAGGAGATTTCTTTGGGAAGT
 TAAAATCACAATAACAAAGGTGCAAAACAACAACACACAGATGATAGTGTTCAGCTCTGCAACAA
 ATACCAACCGAGGCTGCACATCGTGGAAGTGACAGAGGACGGTGTGGAGGACTTGAATGAACCTTCCAAG
 ACTCAGACCTTACCTTCTCAGAGACACAGTTCATCGCTGTGACGGCCTACCAAAACACGGATATCACCC
 AGCTAAAGATCGACCATAACCCCTTCGCCAAAGGCTTCCGGGACAACTACGATTCATGTACACGGCTTC
 AGAAAATGACAGGTTAACTCCATCTCCACGGATTCCCCTAGATCCCATCAGATTGTCCCTGGAGGTCGG
 TACGGCGTTCAAACCTTCTTCCCGGAGCCCTTTGTCAACACTTTGCCTCAAGCCGATATTATAACGGTG
 AGAGAACCGTGCCACAGACCAACGGCTCCTCTACCCCAACAGAGCGAAGAGGTGGCAACCTCCCCA
 GCGGTGGCTTGTACGCCTGTCCAGCAACCTGTGACCAACAAGCTAGACATCGGTTCTTATGAATCTGAA
 TATACTTCCAGTACCTTGTCCCATATGGTATTAAGTCTTGCCCTGCAGACATCCCATGCCCTGGGGT
 ATTACCCTGACCCGACCTTCCCTGCTATGGCAGGGTGGGAGGCCGTGGCGCTTATCAGAGGAAGATGGC
 AGCTGGACTACCATGGACATCCAGAATGAGCCACCTGTCTTCCAGAAGATCAGCTTGCCAAGGAAAAA
 GTTAAAGAAGAGATTAGTTCTCCTGGATAGAGACTCCCCCTCCATCAAGTCTCTAGACTCCAGCGACT
 CCGGGGTGTACAACAGCGCTTGAAGAGAAAGCGCTGTCTCCAGCACCCCCAGCAATGGAAACTCGCC
 CCCATAAAGTGTGAGGACATTAACACTGAAGAGTACAGTAAAGACACCTCCAAGGCATGGGGGCTTAT
 TATGCTTTTTACACAAGTCCC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_010136.3, NP_034266.2</p>
RefSeq Size:	<p>3731 bp</p>
RefSeq ORF:	<p>2124 bp</p>
Locus ID:	<p>13813</p>
UniProt ID:	<p>O54839</p>
Cytogenetics:	<p>9 F3</p>
MW:	<p>75.3 kDa</p>
Gene Summary:	<p>Functions as a transcriptional activator playing a crucial role during development. Functions in trophoblast differentiation and later in gastrulation, regulating both mesoderm delamination and endoderm specification. Plays a role in brain development being required for the specification and the proliferation of the intermediate progenitor cells and their progeny in the cerebral cortex. Also involved in the differentiation of CD8+ T-cells during immune response regulating the expression of lytic effector genes.[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR210179