

Product datasheet for **MR210206L4V**

Aggf1 (NM_025630) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aggf1 (NM_025630) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aggf1
Synonyms:	2010009L17Rik; 2310029P06Rik; AW112072; Peg3; VG5Q
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025630
ORF Size:	2136 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210206).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025630.2
RefSeq Size:	3262 bp
RefSeq ORF:	2136 bp
Locus ID:	66549
UniProt ID:	Q7TN31
Cytogenetics:	13 D1
Gene Summary:	Promotes angiogenesis and the proliferation of endothelial cells. Able to bind to endothelial cells and promote cell proliferation, suggesting that it may act in an autocrine fashion (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)