

Product datasheet for **MR210211L4V**

Ankrd6 (NM_080471) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ankrd6 (NM_080471) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ankrd6
Synonyms:	di
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080471
ORF Size:	2139 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR210211).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080471.3 , NP_536719.2
RefSeq Size:	4573 bp
RefSeq ORF:	2139 bp
Locus ID:	140577
UniProt ID:	Q69ZU8
Cytogenetics:	4 A5



[View online »](#)

Gene Summary:

This gene encodes a protein which is thought to be involved in the Wnt signaling pathway and embryonic axis formation. Similar genes have been found in human, rhesus macaque, and zebrafish. Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]